

Cardinal Health™ COOLSOURCE™ Cooling System

Instructions for Use

The Cardinal Health™ COOLSOURCE™ Cooling Vest is intended for use with Cardinal Health™ Cooling Packs to provide cooling and comfort to healthcare professionals.

The COOLSOURCE™ Cooling System is non-sterile and should be worn underneath a sterile surgical gown or other protective infection control apparel as required for each procedure.

1. Store the COOLSOURCE™ Cooling Vest in a cool, dry place with other non-sterile personal protective equipment.
2. Allow at least two hours for the COOLSOURCE™ Cooling Packs to freeze before use. For best results freeze the COOLSOURCE™ Cooling Packs for at least 24 hours before the first use for maximum efficacy.

Note: The COOLSOURCE™ Cooling Vest should not be placed in the freezer with the Cooling Packs.

3. Insert COOLSOURCE™ Cooling Packs into vest pockets.

Note: Use of COOLSOURCE™ Cooling Packs packs can be customized for each user. If application is not indicated in all six areas (such as the shoulder area), they can be omitted.

4. Don the COOLSOURCE™ Cooling Vest before required protective apparel.

Note: If the COOLSOURCE™ Cooling Vest is in a sterile procedure pack it cannot be donned using sterile technique. Once removed from the sterile pack, the COOLSOURCE™ Cooling Packs should be inserted into vest pockets and donned before scrubbing and donning the sterile surgical gown.

5. Tie belt in the front to optimize fit and positioning of COOLSOURCE™ Cooling Packs.
6. Proceed to don required personal protective equipment.

Removal

1. Remove personal protective equipment.
2. Remove COOLSOURCE™ Cooling Packs from COOLSOURCE™ Cooling Vest.
 - Inspect COOLSOURCE™ Cooling Packs. If they are soiled or worn, dispose of them.
 - If not soiled, Cooling Packs should be wiped with the facility's disinfecting wipes before returning to freezer for re-use.*
 - Follow your facility's infection control guidelines for re-use of supplies.
3. Dispose of COOLSOURCE™ Cooling Vest per facility guidelines.

*The COOLSOURCE™ Cooling Packs can be re-used until soiled or there are visible signs of wear and tear. Follow your facility's infection control protocol for the re-use of OR supplies or the protocol in other departments where COOLSOURCE™ is used.