

 **SAFE, EFFECTIVE WOUND CLEANSING**
 **NON-CYTOTOXIC**
 **NO SODIUM HYPOCHLORITE (BLEACH)**

BIHOCL™ PureCleanse™ is an effective-yet-gentle wound cleansing solution that can be used one or more times each day to help keep your wound and the surrounding area clean and reduce bioburden.

BIHOCL™ PureCleanse™ is a pure, stable, scientifically formulated >99.9% pure hypochlorous acid solution, used to treat acute and chronic wounds, burns and post-surgical wounds.



**A MADE-IN-CANADA SOLUTION
RECOMMENDED BY DOCTORS & NURSES ACROSS CANADA**

 **HIGHLY EFFECTIVE**
 **STING-FREE**
 **EASY TO APPLY**

 **ADDITIVE-FREE**
 **OIL-FREE**
 **ALCOHOL-FREE**

 **VEGAN & ORGANIC**
 **ECO-SAFE**
 **ALL NATURAL**

Patient Care Instructions

Always follow directions provided by your healthcare professional. Below are general guidelines for using the product.

- 1 Apply BIHOCL™ PureCleanse™ solution liberally to the wound and surrounding skin. Can also be used as a soak for 5-10 minutes.**
- 2 When application or soak is complete, pat the area dry with sterile gauze.**
- 3 Re-apply wound dressings as directed by your healthcare professional.**

See related Instructions for Use (IFU) documentation for detailed instructions.

What is Hypochlorous?

Hypochlorous acid (HOCl) is a powerful antimicrobial component of our bodies' internal defense system, produced by our white blood cells. HOCl is the natural molecule contained within BIHOCL™ PureCleanse™ which is safe, effective, and mimics the body's natural response to harmful pathogens.

BIHOCL™ PureCleanse™ is the only hypochlorous formulation tested by the National Institute of Health (NIH) for purity and is unsurpassed in its safety and stability.

- Strictly controlled pH
- Class II Health Canada Medical Device
- 2 year shelf life, 30 days once opened

Ordering Information



PS4018 PureCleanse™ Antimicrobial Cleansing Solution (33055-118DC)
118 mL / EA

PS4019 PureCleanse™ Antimicrobial Cleansing Solution (33055-250DC)
250 mL / EA

BIOMIQ

Biomiq is a Canadian medical device company on a mission to improve patient outcomes through innovations that harness the genius of the natural world around us.

1-800-573-1295 | info@biomiq.health | www.Biomiq.health

Manufactured for Biomiq Inc. by Sterasure Inc.
137 Glasgow St., Unit 115, Kitchener, ON, N2G 4X8, Canada

Rev004