



Chemotherapy Protection and Safety Products Catalogue



Table of Contents

	Protective Gloves 3
	ChemoPlus™ Exam Gloves - Latex
	ChemoPlus™ Neoprene Exam Gloves - Latex Free
	ChemoPlus™ Nitrile Exam Gloves - Latex Free
	Protective Apparel 5
	ChemoPlus™ Poly-Coated Chemotherapy Gowns
	ChemoBloc™ Poly-Coated Chemotherapy Gowns
	Poly-coated SMS Chemotherapy Gowns
	Kendall™ Splash Resistant Gowns
	ChemoPlus™ Protective Sleeves
	Kendall™ Protective Lab Coat
	Kendall™ Smocks
	Dust/Mist Particulate Respirator Mask
	ChemoPlus™ Protective Eyewear
	Kits 9
	ChemoPlus™ Preparation and Administration Kits
	ChemoPlus™ Spill Kits
	BioBloc™ Spill Kit
	Absorbent Products 12
	ChemoPlus™ Absorbent Pads
	ChemoPlus™ Low Lint Towels
	ChemoPlus™ Prep Mats
	ChemoPlus™ Floor Mats
	Waste Containers, Bags and Labels 13
	ChemoPlus™ Chemo Soft Waste Container
	ChemoPlus™ Waste Bags
	ChemoPlus™ Drug Transport Bags
	ChemoPlus™ Medical Labels
	Chemotherapy Sharps Containers 14
	ChemoSafety™ Sharps Container with Slide Lid or Hinged Lid
	IVA Security Seals 15
	ChemoPlus™ IVA Seals for Bottles & Vials
	ChemoPlus™ IVA Seals for I.V. Bags
	ChemoPlus™ IVA Seals for Syringes & Medication Containers
	Training Program 17
	ChemoPlus™ Chemotherapy Training Kit
	Vial Venting System 17
	ChemoPlus™ Vial Vent Chemo Dispensing Pin
	Index 18

Protective Gloves

ChemoPlus™ Exam Gloves - Latex

- 18 mil thickness at fingertips
- Glove length from tip to cuff approximately 30.5 cm (12 in.) to protect wrists
- Beaded cuff resists tearing
- Recommended for maximum protection when preparing chemotherapy drugs or cleaning chemotherapy spills
- Hand specific and textured
- Contains latex
- Powder free



ZCT0194-1

Item Code	Description	Size	Thickness	Packaging
ZCT0191-1	Latex Exam Gloves	Small	18 mil	50 pr/bx, 10 bx/ca
ZCT0192-1	Latex Exam Gloves	Medium	18 mil	50 pr/bx, 10 bx/ca
ZCT0193-1	Latex Exam Gloves	Large	18 mil	50 pr/bx, 10 bx/ca
ZCT0194-1	Latex Exam Gloves	X-Large	18 mil	25 pr/bx, 10 bx/ca

ChemoPlus™ Exam Gloves - Latex

- 10 mil thickness at fingertips
- Glove length from tip to cuff approximately 30.5 cm (12 in.) to protect wrists
- Beaded cuff resists tearing
- Provides maximum protection when administering chemotherapy drugs
- Hand specific and textured
- Contains latex
- Powder free



ZCT5055G

Item Code	Description	Size	Thickness	Packaging
ZCT5055G	Latex Exam Gloves	Small	10 mil	50 pr/bx, 4 bx/ca
ZCT5056G	Latex Exam Gloves	Medium	10 mil	50 pr/bx, 4 bx/ca
ZCT5057G	Latex Exam Gloves	Large	10 mil	50 pr/bx, 4 bx/ca

ChemoPlus™ Neoprene Exam Gloves - Latex Free

- Softer and more comfortable than other latex free gloves
- 9 mil neoprene gloves offer greater protection than latex gloves of comparable thickness
- Glove length from tip to cuff approximately 30.5 cm (12 in.) to protect wrists
- Beaded cuff resists tearing
- Recommended for maximum protection when preparing chemotherapy drugs or cleaning chemotherapy spills
- Hand specific and textured
- Powder free



ZCT21921

Item Code	Description	Size	Thickness	Packaging
ZCT21911	Neoprene Exam Gloves	Small	9 mil	50 pr/bx, 3 bx/ca
ZCT21921	Neoprene Exam Gloves	Medium	9 mil	50 pr/bx, 3 bx/ca
ZCT21931	Neoprene Exam Gloves	Large	9 mil	50 pr/bx, 3 bx/ca
ZCT21941	Neoprene Exam Gloves	X-Large	9 mil	50 pr/bx, 3 bx/ca

ChemoPlus™ Nitrile Exam Gloves (Non-Sterile and Sterile) - Latex Free

- Gloves are ambidextrous
- 8 mil thickness at fingertips
- Unique tear-away feature indicates when glove has been punctured
- Glove length from tip to cuff approximately 30.5 cm (12 in.) to protect wrists
- Beaded cuff resists tearing
- Provides maximum protection when administering chemotherapy drugs
- Textured
- Powder free



ZCT5072G

Item Code	Description	Size	Thickness	Packaging
Non-Sterile				
ZCT5071G	Nitrile Exam Gloves	Small	8 mil	100 ea/bx, 4 bx/ca
ZCT5072G	Nitrile Exam Gloves	Medium	8 mil	100 ea/bx, 4 bx/ca
ZCT5073G	Nitrile Exam Gloves	Large	8 mil	100 ea/bx, 4 bx/ca
ZCT5074G	Nitrile Exam Gloves	X-Large	8 mil	100 ea/bx, 4 bx/ca
Sterile				
ZCTS731S	Nitrile Exam Gloves	Small	8 mil	50 pr/bx, 4 bx/ca
ZCTS731M	Nitrile Exam Gloves	Medium	8 mil	50 pr/bx, 4 bx/ca
ZCTS731L	Nitrile Exam Gloves	Large	8 mil	50 pr/bx, 4 bx/ca
ZCTS731XL	Nitrile Exam Gloves	X-Large	8 mil	50 pr/bx, 4 bx/ca

Sterile gloves are packaged in non-linting peel pouch, 2 per pouch.

Protective Apparel

ChemoPlus™ Poly-Coated Chemotherapy Gowns - Maximum Protection

- Splash resistant poly-coated front and sleeves
- Lightweight, low linting fabric
- Closed front and wrap-around style for added protection against spills
- Closed back design
- Long sleeves with knit cuffs
- Waist ties with Velcro neck closures
- Meets OSHA, ASHP, and ONS recommendations for wearing a gown while preparing or administering chemotherapy drugs
- Latex free

Item Code	Description	Size	Colour	Packaging
ZCT5502	Poly-Coated Chemotherapy Gowns	Medium	Blue	30 ea/ca
ZCT5503	Poly-Coated Chemotherapy Gowns	Large	Blue	10 ea/pk, 3 pk/ca
ZCT5504	Poly-Coated Chemotherapy Gowns	X-Large	Blue	30 ea/ca
ZCT5505	Poly-Coated Chemotherapy Gowns	XX-Large	Blue	30 ea/ca



ZCT5502

ChemoBloc™ Poly-Coated Chemotherapy Gowns - Maximum Protection

- Open back design
- Splash resistant poly-coated front and sleeves
- Lightweight, low linting fabric
- Closed front and side design for added protection against spills
- Long sleeves with knit cuffs
- Waist ties with Velcro neck closures
- Open back for maximum ventilation
- Meets OSHA, ASHP, and ONS recommendations for wearing a gown while preparing or administering chemotherapy drugs
- Latex free



Back view

Item Code	Description	Size	Colour	Packaging
ZDP5003G	Poly-Coated Chemotherapy Gowns	Medium	Light Blue	24 ea/ca
ZDP5001G	Poly-Coated Chemotherapy Gowns	Large	Light Blue	12 ea/bg, 2 bg/ca
ZDP5002G	Poly-Coated Chemotherapy Gowns	X-Large	Light Blue	24 ea/ca
ZDP5004G	Poly-Coated Chemotherapy Gowns	XX-Large	Light Blue	24 ea/ca



ZDP5003G

Poly-coated SMS Chemotherapy Gowns

- Engineered to provide users with protection from the potential hazards associated with handling chemotherapy agents
- Poly-coated material for durability and protection
- Lightweight, closed-back design, secured with ties for maximum comfort
- Low-linting material helps to provide a clean environment
- Long sleeves with knit cuffs provide added comfort

Item Code	Description	Size	Colour	Packaging
A8200CG	Poly-Coated SMS Chemotherapy Gowns	Universal	Blue	10 ea/pk, 10 pk/ca
A8201CG	Poly-Coated SMS Chemotherapy Gowns	X-Large	Blue	10 ea/pk, 10 pk/ca

Kendall™ Splash Resistant Gowns

- Lightweight, low linting
- Closed front and wrap-around style for added protection
- Long sleeves with knit cuffs
- Waist ties with Velcro™ neck closures
- Latex free

Item Code	Description	Size	Colour	Packaging
ZCT5000	Splash Resistant Protective Gown	Medium	Blue	10 ea/pk, 3 pk/ca
ZCT5100	Splash Resistant Protective Gown	Large	Blue	10 ea/pk, 3 pk/ca
ZCT5101	Splash Resistant Protective Gown	X-Large	Blue	30 ea/ca
ZCT5102	Splash Resistant Protective Gown	XX-Large	Blue	30 ea/ca



ZCT5000

ChemoPlus™ Protective Sleeves

- Splash resistant poly-coated sleeves
- Lightweight, low linting fabric
- Elastic ends to maintain secure non-slip fit
- Knit cuffs
- Universal fit
- Latex free

Item Code	Description	Size	Colour	Packaging
ZCT5500	Poly-Coated Sleeves	Universal	Blue	100 ea/bx, 6 bx/ca



ZCT5500

Kendall™ Protective Lab Coat

- General use lab coat
- Long sleeves with knit cuffs
- Convenient front pocket
- Lapel can be worn open or closed
- Latex free

Item Code	Description	Size	Colour	Packaging
ZBB6502L	Protective Lab Coat	Medium	White	20 ea/ca
ZBB6503L	Protective Lab Coat	Large	White	20 ea/ca
ZBB6504L	Protective Lab Coat	X-Large	White	20 ea/ca



ZBB6502L

Kendall™ Smocks

- General use smocks
- Long sleeves with knit cuffs
- Convenient front pocket
- Lapel can be worn open or closed

Item Code	Description	Size	Colour	Packaging
Z932879	Smock with Collar	Small	Teal	24 ea/ca
Z932882	Smock with Collar	Medium	Teal	24 ea/ca
Z932885	Smock with Collar	Large	Teal	24 ea/ca
Z932886	Smock with Collar	X-Large	Teal	24 ea/ca



Z932879

Dust/Mist Particulate Respirator Mask

- N95 Rated
- Resists collapsing
- Easy donning
- Pre-fitted nose bridge
- Tested and certified by NIOSH to have a filter efficiency level of 95% or greater against solid and non-liquid (non-oil based) aerosols
- Meets CDC guidelines for TB exposure control standards and protection against laser plumes or surgical smoke
- Bacterial filtration efficiency (BFE) > 99.9%
- Latex free



ZCT0600-1

Item Code	Description	Size	Colour	Packaging
ZCT0600-1	Dust/Mist Particulate Respirator Mask	Universal	White	20 ea/bx, 6 bx/ca

ChemoPlus™ Protective Eyewear

- Splash resistant
- Universal fit
- Comfortably fits over most eyeglasses
- Ventilated



ZDP5030G

Item Code	Description	Size	Packaging
ZDP5030G	Protective Wrap-Around Goggles	Universal	6 ea/ca
ZCT0400-1	Plastic Eyewear with Side Shields	Universal	10 ea/bx, 5 bx/ca



ZCT0400-1

Kits

ChemoPlus™ Preparation and Administration Kits

Item Code	Description	Packaging
ZCT4010	Kit Includes: Gown - 1 ChemoPlus™ Poly-Coated Chemotherapy Gown, Large Latex Exam Glove - 1 pair Medium, 10 mil, Textured Transport Bag - 1 SafeLock™ Bag (12 x15 in.) Signs - 2 "Caution" and "Flush Twice" Prep Mats - 1 (11 x 17 in.)	24 ea/ca
ZCT4020	Kit Includes: Gown - 1 ChemoPlus™ Poly-Coated Chemotherapy Gown, Large Latex Exam Glove - 1 pair Medium, 18 mil, textured Transport Bag - 1 SafeLock™ Bag (12 x 15 in.) Signs - 2 "Caution" and "Flush Twice" Prep Mats - 1(11 x 17 in.)	24 ea/ca
ZCT4025	Kit Includes: Gown - 1 ChemoPlus™ Poly-Coated Chemotherapy Gown, X-Large Latex Exam Glove - 1 pair Medium, 10 mil, Textured Transport Bag - 1 SafeLock™ Bag (12 x15 in.) Signs - 2 "Caution" and "Flush Twice" Prep Mats - 1 (11 x17 in.)	24 ea/ca
ZCT4100	Kit Includes: Gown - 1 ChemoPlus™ Poly-Coated Chemotherapy Gown, Medium Latex Exam Glove - 1 pair Small & 1 pair Medium, 18 mil, Textured Mask - 1 Dust/Mist Respirator Mask Transport Bag - 1 SafeLock™ Bag (12 x 15 in.) Prep Mats - 1 (11 x 17 in.)	24 ea/ca
ZCT4101	Kit Includes: Gown - 1 ChemoPlus™ Poly-Coated Chemotherapy Gown, Large Latex Exam Glove - 1 pair Medium & 1 pair Large, 18 mil, Textured Mask - 1 Dust/Mist Respirator Mask Transport Bag - 1 SafeLock™ Bag (12 x 15 in.) Prep Mats - 1 (11 x 17 in.)	24 ea/ca

ChemoPlus™ Spill Kits

- Designed for safe and fast clean-up of chemotherapy drug spills
- Packaged to facilitate compliance with OSHA and ASHP guidelines

Item Code	Description	Packaging
ZDP5108K	Kit Includes: Gown - 1 Poly-Coated Chemotherapy Gown, Maximum Protection, Blue, XX-Large Latex Free Nitrile Gloves - 2 pair Medium, 8 mil Mask - 1 Dust/Mist Respirator Mask Waste Bag - 1 (15 gal.) Chemo Bio-Wipe™ Bag - 1 Goggles - 1	24 ea/ca
ZDP5016K	Kit Includes: Gown - 1 Poly-Coated Chemotherapy Gown, Maximum Protection, Blue, X-Large Latex Free Nitrile Gloves - 1 pair Medium & 1 pair Large, 8 mil Mask - 1 Dust/Mist Respirator Waste Bag - 1 (2 gal.), 1 (15 gal.) Shoe Covers - 1 pair Goggles - 1 pair Wrap-around Goggles Signs - 1 "Caution: Dry Spill" Scoop & Brush - 1 set Labels - 6 "Caution Chemotherapy Waste" and "Drug Spill" Absorbent Towels - 3 (12 x12 in.) Spill Control Pillows - 2 Report Forms - 1 Exposure Report Form	6 ea/ca
ZCT4004	Kit Includes: Gown - 1 Poly-Coated Chemotherapy Gown, Maximum Protection, Blue, Large Nitrile Gloves - 1 pair Large & 1 pair X-Large Mask - 1 Dust/Mist Respirator Waste Bag - 2 (20 gal.) Goggles - 1 pair Wrap-around Goggles Spill Towel - 3 Shoe Covers - 1 pair Signs - 1 Caution - "Caution: Dry Spill" ChemoSorb™ Pads - 3 Scoop & Scraper - 1 set Tie Wrap - 2 Hazardous Drug Exposure Form - 1	4 ea/ca

BioBloc™ Spill Kit

- Designed to facilitate fast and safe clean-up of blood and body fluid spills
- Can be used in numerous occupations
- Kits comply with OSHA's rule for Occupational Exposure to Bloodborne Pathogens
- Latex free



Item Code	Description	Packaging
ZBB6016K	Kit Includes: Gown - 1 Poly-Coated Chemotherapy Gown, Maximum Protection, X-Large Latex Free Gloves - 1 pair Medium & 1 pair Large, 8 mil Mask - 1 with eye shield & ear loops Bio-Hazard Waste Bag - 1 (10 gal.) & 1 (33 gal.) bags Wipes - 2 Antimicrobial Disposable Cloths - 2 Germicidal Sanizide Plus Packet - 1 Fluid Control Packets - 1 Scoop & Scraper - 1 set Labels - 2 "BioHazard: Handle With Caution" Forms - 1 Incident Report	6 ea/ca

Absorbent Products

ChemoPlus™ Absorbent Pads

- Highly absorbent
- ZBS0013 pads are impregnated with Green-Z™ solidifying agent to instantly turn liquid spills into a semi-gel
- Latex free, single use


Item Code	Description	Size	Packaging
ZBS0013	Absorbent Pad	12 x 12 in. (30.5 x 30.5 cm)	100 ea/ca



ZBS0013

ChemoPlus™ Low Lint Towels

Item Code	Description	Size	Packaging
ZCT0014	Low Lint Towel	9 x 9 in. (23 x 23 cm)	300 ea/ca



ZCT0014

ChemoPlus™ Prep Mats

- Highly absorbent design
- Available in several sizes
- Economy version available
- Latex free


Item Code	Description	Size	Packaging
ZDP5010M	Economy Absorbent Prep Mat	17 x 19 in. (43.2 x 48.3 cm)	125 ea/ca
ZCT0300-1	Prep Mat Spill Absorbent	11 x 17 in. (27.9 x 43.2 cm)	50 ea/ca, 5 bx/ca
ZCT0302-1	Prep Mat Spill Absorbent	16 x 22 in. (40.6 x 55.9 cm)	50 ea/ca, 5 bx/ca
ZCT0303-1	Prep Mat Spill Absorbent	16 x 69 in. (40.6 cm x 1.8 m)	15 ea/ca, 5 bx/ca



ZCT0300-1

ChemoPlus™ Floor Mats

- Convenient and economical way to trap impurities from foot traffic in the pharmacy
- Easy to adhere to floor
- Simple to peel back, exposing new layer
- Thirty peel-back layers per mat
- Latex free



ZCT0007

Item Code	Description	Size	Colour	Packaging
ZCT0007	Adhesive Floor Mat	18 x 46 in. (45.7 cm x 1.2 m)	White	4 ea/ca
ZCT0071	Adhesive Floor Mat	18 x 46 in. (45.7 cm x 1.2 m)	Blue	4 ea/ca

Waste Containers, Bags and Labels

ChemoPlus™ Chemo Soft Waste Container

- Easy to assemble
- Dual flap opening with self-locking lid
- Holds approximately 20 gal. of waste
- ZCT1500 includes one 20 gal. waste bag
- Latex free

Item Code	Description	Size	Thickness	Colour	Packaging
ZCT1500	Chemo Soft Waste Container	20 gal.	4 mil	Yellow	6 ea/ca



ZCT1500

ChemoPlus™ Waste Bags

- Available in a variety of sizes
- 4 mil thick, very durable
- Chemotherapy yellow color
- Closure ties included
- Latex free

Item Code	Description	Size	Thickness	Colour	Packaging
ZDP5041B	Chemo Soft Waste Bag	2 gal.	4 mil	Yellow	250 ea/ca
ZDP5042B	Chemo Soft Waste Bag	15 gal.	4 mil	Yellow	100 ea/ca
ZCT2100	Chemo Soft Waste Bag	20 gal.	4 mil	Yellow	100 ea/ca
ZDP5043B	Chemo Soft Waste Bag	30 gal.	4 mil	Yellow	100 ea/ca



ZCT2100

ChemoPlus™ Drug Transport Bags

- Available in two sizes
- 4 mil thick for maximum protection
- Zip-lock closure
- Clear to visualize contents
- Latex free

Item Code	Description	Size	Thickness	Colour	Packaging
ZCT0575	Chemo Drug Transport Bag	6 x 9 in. (15.2 x 22.9 cm)	4 mil	White	200 ea/ca
ZCT0500	Chemo Drug Transport Bag	12 x 15 in. (30.5 x 38.1 cm)	4 mil	White	250 ea/ca



ZCT0575

ChemoPlus™ Medical Labels

- Includes labels for pharmacy, oncology, I.V. therapy, and patient care
- Reduces risk of medication errors and ensures clinician and patient safety
- Pre-printed message provides convenient and cost-effective communication
- "Caution: Chemotherapy Handle with gloves, Dispose of properly"
- Latex free

Item Code	Description	Size	Packaging
ZCT6006-1	Chemo Caution Label	1/2 x 3 in. (1.3 x 7.6 cm)	500 ea/bx, 4 bx/ca
ZCT6011-1	Chemo Caution Label	1 x 3 in. (2.5 x 7.6 cm)	500 ea/bx, 4 bx/ca



ZCT6011-1

Chemotherapy Sharps Containers

ChemoSafety™ Sharps Containers with Slide Lid or Hinged Lid

- For the disposal of trace chemotherapy waste
- Sealing gasket lids provide leak-resistant containment of fluid
- Absorbent pad in bottom of container absorbs excess liquid
- Adjustable sliding lid reduces the possibility of aerosolization and accommodates chemotherapy waste of various types and sizes
- Available carts with wheels and built-in carrying handles for easy transport and lifting of heavy, full containers
- Containers provide temporary and final locking closure for disposal
- Nestable containers save valuable storage space
- Sliding gasket lids have a round 6 in. opening to accommodate up to 1000 mL bottles
- Bilingual label (French, English)



Z8934CRS



Z8939CRS



Z8985CRS



Z8982CRH



Z8985CRH

Item Code	Description	Size (H x D x W)	Colour	Packaging
Sliding Lid				
Z8985CRS	Sharps Container	8 gal. (17.75 x 11 x 15.5 in.)	Red	10 ea/ca
Z8934CRS	Sharps Container	12 gal. (18.75 x 12.75 x 18.5 in.)	Red	10 ea/ca
Z8939CRS	Sharps Container	18 gal. (26 x 12.75 x 18.25 in.)	Red	5 ea/ca
Hinge Lid				
Z8982CRH	Sharps Container	2 gal. (10 x 7.25 x 10.5 in.)	Red	20 ea/ca
Z8985CRH	Sharps Container	8 gal. (17.75 x 11 x 15.5 in.)	Red	10 ea/ca
Z8989CRH	Sharps Container	18 gal. (26 x 12.75 x 18.25 in.)	Red	5 ea/ca

IVA Security Seals

ChemoPlus™ IVA Seals for Bottles & Vials

- Dual warning system: Seals self-destruct upon attempted removal and needle punctures are highly visible in laminated foil
- Sterile IVA seals protect point of entry of IV admixtures from contamination and accidental double dosing
- Multiple sizes and colours available to fit most bottles and vials
- Latex free & PVC free



Item Code	Description	Packaging
For Top Vials		
ZCP3005A	Silver IVA Seals, 9 mm	1100 ea/bx, 12 bx/ca
ZCP3005R	Red IVA Seals, 9 mm	1100 ea/bx, 12 bx/ca
ZCP3003A	Silver IVA Seals, 13 mm	1100 ea/bx, 20 bx/ca
ZCP3003B	Blue IVA Seals, 13 mm	1100 ea/bx, 20 bx/ca
ZCP3003G	Green IVA Seals, 13 mm	1100 ea/bx, 20 bx/ca
ZCP3003R	Red IVA Seals, 13 mm	1100 ea/bx, 20 bx/ca
ZCP3003Y	Yellow IVA Seals, 13 mm	1100 ea/bx, 20 bx/ca
ZCP3002G	Green IVA Seals, 20 mm	1000 ea/bx, 12 bx/ca
ZCP3002B	Blue IVA Seals, 20 mm	1000 ea/bx, 12 bx/ca
ZCP3002Y	Yellow IVA Seals, 20 mm	1000 ea/bx, 12 bx/ca
ZCP3002O	Orange IVA Seals, 20 mm	1000 ea/bx, 12 bx/ca
ZCP3002R	Red IVA Seals, 20 mm	1000 ea/bx, 12 bx/ca
For Top Bottles & Piggyback Containers		
ZCP3001A	Silver IVA Seals, 28 mm	1000 ea/bx, 12 bx/ca
ZCP3001G	Green IVA Seals, 28 mm	1000 ea/bx, 12 bx/ca
ZCP3001B	Blue IVA Seals, 28 mm	1000 ea/bx, 12 bx/ca
ZCP3001R	Red IVA Seals, 28 mm	1000 ea/bx, 12 bx/ca
For Top Bottles & Lily's Faspak		
ZCP3009A	Silver IVA Seals, 36 mm	1000 ea/bx, 12 bx/ca
ZCP3009G	Green IVA Seals, 36 mm	1000 ea/bx, 12 bx/ca
ZCP3009R	Red IVA Seals, 36 mm	1000 ea/bx, 12 bx/ca


Most IVA seals are available in these six colors:



A-Silver



R-Red




B-Blue



G-Green



O-Orange



Y-Yellow

ChemoPlus™ IVA Seals for I.V. Bags

- Visual communication that prevents dosing errors
- Dual warning system: Seals self-destruct upon attempted removal and needle punctures are highly visible in laminated foil
- IVA seals on bags instantly alert staff that a drug has been added
- Sterile IVA seals protect point of entry of IV admixtures from contamination and accidental double dosing
- Three styles to accommodate most commonly used IV bag
- PVC free
- Latex free



Item Code	Description	Packaging
For Hospira's Large Bags & Add-Vantage™ System		
ZCP3004B	Blue IVA Seals	1000 ea/bx, 12 bx/ca
For Baxter's Viaflex™ & Mini-Bag Plus Containers		
ZCP3011A	Silver IVA Seals	1000 ea/bx, 20 bx/ca
ZCP3011G	Green IVA Seals	1000 ea/bx, 20 bx/ca
ZCP3011B	Blue IVA Seals	1000 ea/bx, 20 bx/ca
ZCP3011R	Red IVA Seals	1000 ea/bx, 20 bx/ca
ZCP3011Y	Yellow IVA Seals	1000 ea/bx, 20 bx/ca
For B. Braun's Excel™ Bags & Abbott's Small Bags		
ZCP3013G	Green IVA Seals	1000 ea/bx, 12 bx/ca
ZCP3013B	Blue IVA Seals	1000 ea/bx, 12 bx/ca
ZCP3013R	Red IVA Seals	1000 ea/bx, 12 bx/ca



ZCT3004B



ZCT3013R

ChemoPlus™ IVA Seals for Syringes & Medication Containers

- Ensures the security of needle covers, syringe caps, box tops and many other medication containers
- Non-sterile, tamper evident seals self-destruct upon attempted removal
- Unused medications may be returned and reissued with confidence if applied seals are intact
- Available in large, 4-1/4 in. (10.8 cm) and small, 3-5/8 in. (9.2 cm) lengths



Item Code	Description	Packaging
For Larger Syringes, 4-1/4 in. length		
ZCP3020B	Blue IVA Seals	1000 ea/bx, 12 bx/ca
ZCP3020Y	Yellow IVA Seals	1000 ea/bx, 12 bx/ca
ZCP3020R	Red IVA Seals	1000 ea/bx, 12 bx/ca
For Smaller Syringes, 3-5/8 in. length		
ZCP3030R	Red IVA Seals	1000 ea/bx, 12 bx/ca



ZCP3020B

Training Program

For Pharmacy & Oncology Healthcare Providers

The ChemoChek Training and Certification program provides participants with self-administered, hands-on training for the techniques required to safely handle chemotherapy drugs. At the conclusion of this program, the participant will be able to understand and use proper procedures while handling cytotoxic drugs.

ChemoPlus™ Chemotherapy Training Kit

- Educates participants about the risks associated with the handling of chemotherapy drugs
- Refill kits available
- Note: One certificate of completion will be issued for each kit purchased
- Latex free



Item Code	Description	Packaging
ZCT4400	Pharmacy Training Kit includes: Fluorescein Dye - 1 Set Vials/Ampoule UV Light without Batteries - 1 Training Manual with Certificate Application - 1 Pharmacy Manual Prep Mat - 5 (10 x 17 in.) Transport Bags - 5	1 ea/ca
ZCT4220	Fluorescein Refill Kit includes: Fluorescein Dye - 1 Set Vials/Ampoule	4 ea/ca

ZCT4400



ZCT4220

Vial Venting System

A system designed to prevent aerosolization of chemotherapy and other hazardous drugs into the atmosphere during reconstitution.

ChemoPlus™ Vial Vent Chemo Dispensing Pin

- Can be used for reconstituting chemotherapy drugs, antibiotics, and other hazardous drugs
- 0.22 micron hydrophobic filter to safely remove any aerosolized drugs
- Safe, needleless system prevents needle stick injury
- Unique cannula is easy to insert
- Facilitates rapid transfer of solution
- Sterile, non-pyrogenic fluid path
- Luer-lock connection assures secure, safe connection to syringe
- Single use only
- Latex free



Item Code	Description	Packaging
ZDP5201V	Vial Vent Chemo Dispensing Pin	50 ea/bx, 4 bx/ca

Index

Item Code	Product	Page	Item Code	Product	Page
A8200CG	Poly-Coated SMS Chemotherapy Gowns, Universal	6	ZCT0007	Adhesive Floor Mat, White	12
A8201CG	Poly-Coated SMS Chemotherapy Gowns, X-Large	6	ZCT0014	Low Lint Towel	12
ZBB6016K	BioBloc™ Spill Kit	11	ZCT0071	Adhesive Floor Mat, Blue	12
ZBB6502L	Protective Lab Coat, Medium	7	ZCT0500	Chemo Drug Transport Bag, 12 x 15 in.	13
ZBB6503L	Protective Lab Coat, Large	7	ZCT0575	Chemo Drug Transport Bag, 6 x 9 in.	13
ZBB6504L	Protective Lab Coat, X-Large	7	ZCT1500	Chemo Soft Waste Container, 20 gal.	13
ZBS0013	Absorbent Pad	12	ZCT2100	Chemo Soft Waste Bag, 20 gal.	13
ZCP3001A	Silver IVA Seals, 28 mm	15	ZCT4004	ChemoPlus™ Spill Kit	10
ZCP3001B	Blue IVA Seals, 28 mm	15	ZCT4010	ChemoPlus™ Preparation and Administration Kit	9
ZCP3001G	Green IVA Seals, 28 mm	15	ZCT4020	ChemoPlus™ Preparation and Administration Kit	9
ZCP3001R	Red IVA Seals, 28 mm	15	ZCT4025	ChemoPlus™ Preparation and Administration Kit	9
ZCP3002B	Blue IVA Seals, 20 mm	15	ZCT4100	ChemoPlus™ Preparation and Administration Kit	9
ZCP3002G	Green IVA Seals, 20 mm	15	ZCT4101	ChemoPlus™ Preparation and Administration Kit	9
ZCP3002O	Orange IVA Seals, 20 mm	15	ZCT4220	Fluorescein Refill Kit	17
ZCP3002R	Red IVA Seals, 20 mm	15	ZCT4400	Pharmacy Training Kit	17
ZCP3002Y	Yellow IVA Seals, 20 mm	15	ZCT5000	Splash Resistant Protective Gown, Medium	6
ZCP3003A	Silver IVA Seals, 13 mm	15	ZCT5055G	Latex Exam Gloves, Small	3
ZCP3003B	Blue IVA Seals, 13 mm	15	ZCT5056G	Latex Exam Gloves, Medium	3
ZCP3003G	Green IVA Seals, 13 mm	15	ZCT5057G	Latex Exam Gloves, Large	3
ZCP3003R	Red IVA Seals, 13 mm	15	ZCT5071G	Nitrile Exam Gloves, Non-Sterile, Small	4
ZCP3003Y	Yellow IVA Seals, 13 mm	15	ZCT5072G	Nitrile Exam Gloves, Non-Sterile, Medium	4
ZCP3004B	Blue IVA Seals	16	ZCT5073G	Nitrile Exam Gloves, Non-Sterile, Large	4
ZCP3005A	Silver IVA Seals, 9 mm	15	ZCT5074G	Nitrile Exam Gloves, Non-Sterile, X-Large	4
ZCP3005R	Red IVA Seals, 9 mm	15	ZCT5100	Splash Resistant Protective Gown, Large	6
ZCP3009A	Silver IVA Seals, 36 mm	15	ZCT5101	Splash Resistant Protective Gown, X-Large	6
ZCP3009G	Green IVA Seals, 36 mm	15	ZCT5102	Splash Resistant Protective Gown, XX-Large	6
ZCP3009R	Red IVA Seals, 36 mm	15	ZCT5500	Poly-Coated Sleeves, Blue	6
ZCP3011A	Silver IVA Seals	16	ZCT5502	Poly-Coated Chemotherapy Gowns, Medium, Blue	5
ZCP3011B	Blue IVA Seals	16	ZCT5503	Poly-Coated Chemotherapy Gowns, Large, Blue	5
ZCP3011G	Green IVA Seals	16	ZCT5504	Poly-Coated Chemotherapy Gowns, X-Large, Blue	5
ZCP3011R	Red IVA Seals	16	ZCT5505	Poly-Coated Chemotherapy Gowns, XX-Large, Blue	5
ZCP3011Y	Yellow IVA Seals	16	ZCT0191-1	Latex Exam Gloves, Small	3
ZCP3013B	Blue IVA Seals	16	ZCT0192-1	Latex Exam Gloves, Medium	3
ZCP3013G	Green IVA Seals	16	ZCT0193-1	Latex Exam Gloves, Large	3
ZCP3013R	Red IVA Seals	16	ZCT0194-1	Latex Exam Gloves, X-Large	3
ZCP3020B	Blue IVA Seals	16	ZCT0300-1	Prep Mat Spill Absorbent, 11 x 17 in.	12
ZCP3020R	Red IVA Seals	16	ZCT0302-1	Prep Mat Spill Absorbent, 16 x 22 in.	12
ZCP3020Y	Yellow IVA Seals	16	ZCT0303-1	Prep Mat Spill Absorbent, 16 x 69 in.	12
ZCP3030R	Red IVA Seals	16	ZCT0400-1	Plastic Eyewear with Side Shields	8

Item Code	Product	Page
ZCT0600-1	Dust/Mist Particulate Respirator Mask	8
ZCT21911	Neoprene Exam Gloves, Small	4
ZCT21921	Neoprene Exam Gloves, Medium	4
ZCT21931	Neoprene Exam Gloves, Large	4
ZCT21941	Neoprene Exam Gloves, X-Large	4
ZCT6006-1	Chemo Caution Label, 1/2 x 3 in.	13
ZCT6011-1	Chemo Caution Label, 1 x 3 in.	13
ZDP5001G	Poly-Coated Chemotherapy Gowns, Large, Light Blue	5
ZDP5002G	Poly-Coated Chemotherapy Gowns, X-Large, Light Blue	5
ZDP5003G	Poly-Coated Chemotherapy Gowns, Medium, Light Blue	5
ZDP5004G	Poly-Coated Chemotherapy Gowns, XX-Large, Light Blue	5
ZDP5010M	Economy Absorbent Prep Mat	12
ZDP5016K	ChemoPlus™ Spill Kit	10
ZDP5030G	Protective Wrap-Around Goggles	8
ZDP5041B	Chemo Soft Waste Bag, 2 gal.	13
ZDP5042B	Chemo Soft Waste Bag, 15 gal.	13
ZDP5043B	Chemo Soft Waste Bag, 30 gal.	13
ZDP5108K	ChemoPlus™ Spill Kit	10
ZDP5201V	Vial Vent Chemo Dispensing Pin	17
ZCTS731S	Nitrile Exam Gloves, Sterile, Small	4
ZCTS731M	Nitrile Exam Gloves, Sterile, Medium	4
ZCTS731L	Nitrile Exam Gloves, Sterile, Large	4
ZCTS731XL	Nitrile Exam Gloves, X-Large	4
Z8934CRS	Sharps Container, Sliding Lid, 12 gal.	14
Z8939CRS	Sharps Container, Sliding Lid, 18 gal.	14
Z8982CRH	Sharps Container, Hinge Lid, 2 gal.	14
Z8985CRS	Sharps Container, Sliding Lid, 8 gal.	14
Z8985CRH	Sharps Container, Hinge Lid, 8 gal.	14
Z8989CRH	Sharps Container, Hinge Lid, 18 gal.	14
Z932879	Smock with Collar, Small	7
Z932882	Smock with Collar, Medium	7
Z932885	Smock with Collar, Large	7
Z932886	Smock with Collar, X-Large	7

For more information or to order, contact your local
Cardinal Health Canada **sales representative today or call:**
Acute - 1.888.291.5033
Continuing Care - 1.877.878.7778

