



Cardinal Health™ OR Accessories



Efficiency and adaptability are critical in the OR.

We provide a full range of sterile and non-sterile OR accessories such as skin markers, light handle covers and needle counters, offering you quality, performance and value every time. Products can be added to custom packs or can be ordered individually.

Take any guesswork out of OR operations

OR Safety

Ease adherence to safety standards with our drapes, needle counters and other safety products.

OR Essentials

The last thing you want to think about in the operating room is the efficacy of your equipment. Our comprehensive list of must-haves will put your focus on the task at hand, not the quality of the tools you're using. From marking pens to bandages, we have you covered.

OR Lap Sponges and Towels

Our low lint towels generate 99% less lint than standard cotton towels. Reduce cleanup and protect patient health with sponges and towels that are soft, absorbent and x-ray detectable.




Table of Contents

OR Safety 4 – 7

Magnetic Drapes

Devon™ Magnetic Drapes with Neutral Zones	5
Devon™ Reusable Magnetic Drapes	5
Hands Free Transfer Magnetic Drapes	5

Magnetic Trays and Basins

Needle Counters

Devon™ Foam Strip Needle Counters	6
Devon™ Foam Strip and Magnet Needle Counters	6
Devon™ Single Magnet	6
Foam Block Needle Counters	6
Foam Block and Magnet Needle Counters	6
Magnetic Needle Counters - Double Magnet	7

Safety Ancillary

Devon™ Magnetic Needle Finder	7
-------------------------------	---

Sponge Counting Systems

Sponge Counter Bags	7
Sponge Counting Racks and Brackets	7

Transfer Products

Devon™ Blade Shield Scalpel Holder	7
------------------------------------	---

OR Essentials 8 – 16

Anti-Fog Solution

Devon Surgical Clamps

Disposable Drape Clamps	9
-------------------------	---

Cautery Tip Cleaning Systems

Esmark Bandages

Insufflation Tubing

Light Handle Systems

Universal Light Handles (ULH) and Universal Light Adapters (ULA)	11
Flexible Light Handle Covers	11
Rigid Light Handle	11
Semi-Rigid Light Sleeve	11

Marking Pens, Surgical

Mini Marking Pens	12
Prep Resistant Surgical Marking Pens	12
Traditional Surgical Marking Pens	13

Marking Pens, Surgical Dual Tip

Marking Pens, Utility

Medication Labels

OR Plastics

Surgical Blades

Disposable Surgical Blades	16
----------------------------	----

Surgical Vessel Loops

Mini Surgical Vessel Loops	16
Maxi Surgical Vessel Loops	16

Miscellaneous

OR Lap Sponges and Towels 17 – 19

Lap Sponges

Curity™ Lap Sponges	18
Kerlix™ Lap Sponges	18

OR Towels


Standard OR Towels	19
X-ray Detectable OR Towels	19

Universal Low Lint Towels



For more information or to order, contact your local Cardinal Heath Canada sales representative today or call 1.888.291.5033

OR Safety





Z31140612

Magnetic Drapes

Devon™ Magnetic Drapes with Neutral Zones

- Disposable neutral zone magnetic drapes have a convenient, built in, non-magnetic area to hold or safely transfer sharp instruments
- Sterile and non-sterile options available
- Not made with natural rubber latex

Item No.	Description	Packaging
Z31140612	10 x 16 in. with neutral zone - Devon™	30 ea/ca
Z31140539	16 x 20 in. with neutral zone - Devon™	20 ea/ca
Z31140430	16 x 20 in. with two neutral zones - Devon™	16 ea/ca



Z31140588

Devon™ Reusable Magnetic Drapes

- Constructed from high quality medical grade siliconized rubber for strength and durability
- Round ceramic magnets help prevent perforations after repeated handling and use
- May be ETO or steam-sterilized

Item No.	Description	Packaging
Z31140588	10 x 16 in. reusable - Devon™	3 ea/ca, 1 ca/bx

Disposable Magnetic Drapes

- Disposable magnetic drapes designed to attract and hold surgical instruments on the sterile field
- Helps reduce the risk of percutaneous injuries and helps enhance safe, effective instrument transfer in the OR
- Flexible drapes conform to any shape with non-slip foam backing for secure positioning
- Not made with natural rubber latex

Item No.	Description	Packaging
100-10B	3 x 4 in.	200 ea/ca
200-16B	10 x 16 in.	48 ea/ca
Z31140570	10 x 16 in. - Devon™	30 ea/ca
200-20B	16 x 20 in.	24 ea/ca
Z31140562	16 x 20 in. - Devon™	20 ea/ca



Magnetic Trays and Basins

Devon™ Magnetic Trays and Basins

- Kidney shape basin and rectangular tray easily attach to magnetic drapes
- Not made with natural rubber latex

Item No.	Description	Packaging
Z31140398	3.75 x 9 x 1 in. tray with adhesive strip	50 ea/ca

Needle Counters

- Available in foam block, foam strip, and magnet combinations with needle counts from 10 - 100
- BoxLocks feature provides secure closure and disposal, but may be reopened as needed for recounts
- Blade disarming system simply and safely removes blades for proper disposal
- Adhesive tab allows needle counter to be secured to any work surface
- Not made with natural rubber latex

Devon™ Foam Strip Needle Counters



Z31142170

Item No.	Description	Packaging
With BoxLocks		
Z31142212	15 count/15 capacity, clear top - Devon™	12 ea/bx, 8 bx/ca
Z31142303	30 count/ 30 capacity- Devon™	12 ea/bx, 8 bx/ca
Z31142204	60 count/60 capacity - Devon™	8 ea/bx, 8 bx/ca

Devon™ Foam Strip and Magnet Needle Counters



Z31142410

Item No.	Description	Packaging
With BoxLocks		
Z31142170	15 count/30 capacity - Devon™	12 ea/bx, 8 bx/ca
With Boxlocks and blade disarmer		
Z31142295	15 count/30 capacity, foam strip and magnet, with BoxLock and blade disarmer - Devon™	12 ea/bx, 8 bx/ca

Devon™ Single Magnet

Item No.	Description	Packaging
Z31142410	10 count/10 capacity, clear top - Devon™	21 ea/bx, 8 bx/ca
With BoxLocks		
Z31142139	30 count/30 capacity, clear top - Devon™	8 ea/bx, 8 bx/ca

Foam Block Needle Counters



Z31142469

Item No.	Description	Packaging
3FF80SSA	40 count/80 capacity	6 ea/bx, 8 bx/ca
With BoxLocks		
Z31142485	10 count/20 capacity, adhesive backing - Devon™	12 ea/bx, 8 bx/ca
Z31142501	20 count/40 capacity, adhesive backing - Devon™	12 ea/bx, 8 bx/ca
Z31142402	40 count/70 capacity, adhesive backing - Devon™	8 ea/bx, 8 bx/ca
Z31142469	80 count/80 capacity, double foam - Devon™	64 ea/ca
With BoxLocks and blade disarmer		
Z31142360	20 count/40 capacity, adhesive backing - Devon™	12 ea/bx, 8 bx/ca
Z31142451	30 count/30 capacity - Devon™	12 ea/bx, 8 bx/ca

Foam Block and Magnet Needle Counters

Item No.	Description	Packaging
3FM20SSA	10 count/20 capacity	7 ea/bx, 8 bx/ca
With BoxLocks		
Z31142428	10 count/20 capacity - Devon™	12 ea/bx, 8 bx/ca
Z31142493	20 count/40 capacity - Devon™	12 ea/bx, 8 bx/ca
Z31142311	40 count/70 capacity - Devon™	8 ea/bx, 8 bx/ca
With BoxLocks and blade disarmer		
Z31142394	10 count/20 capacity - Devon™	12 ea/bx, 8 bx/ca
Z31142386	20 count/40 capacity - Devon™	12 ea/bx, 8 bx/ca
Z31142378	50 count/50 capacity - Devon™	12 ea/bx, 8 bx/ca
4FM100SBBL	60 count/100 capacity	40 ea/ca
Z31142352	100 count/100 capacity - Devon™	8 ea/bx, 8 bx/ca



Z31140570

Magnetic Needle Counters - Double Magnet

Item No.	Description	Packaging
3MM20SSA	10 count/20 capacity	12 ea/bx, 6 bx/ca
3MM30SSA	15 count/30 capacity	12 ea/bx, 6 bx/ca
Z31142196	20 count/ 20 capacity - Devon™	21 ea/bx, 8 bx/ca
3MM60SSA	30 count/60 capacity	12 ea/bx, 6 bx/ca
With BoxLocks		
Z31142147	30 count/30 capacity - Devon™	12 ea/bx, 8 bx/ca
Z31142154	60 count/60 capacity - Devon™	8 ea/bx, 8 bx/ca
With BoxLocks and blade disarmer		
Z31142113	30 count/ 30 capacity - Devon™	12 ea/bx, 8 bx/ca



Safety Ancillary

Devon™ Magnetic Needle Finder

- Reusable, rolling magnetic needle finder provides a convenient, safe, and effective method to find and retrieve fallen needles and blades
- Convenient handle length of 44 inches, for ease of use
- Magnets not sold separately
- Non-sterile
- Not made with natural rubber latex

Item No.	Description	Packaging
Z31140422	Rolling Magnetic Needle Finder - Devon™	1 ea/ca

Sponge Counting Systems

- Economical method to develop a standard sponge counting protocol
- Provides a reliable system to assist in the prevention of retained sponges
- Transparent polyethylene pockets with varying backing options provide clear visibility of used sponges
- Racks for I.V. poles include mounting bracket
- Not made with natural rubber latex

Sponge Counter Bags

Item No.	Description	Packaging
Z31144127	10 pocket sponge counter bag dispenser box, white back - Devon™	50/bx, 5 bx/ca
Z5096	10 pocket sponge counter bag dispenser box, blue back - Devon™	50/bx, 5 bx/ca

Sponge Counting Racks and Brackets

Item No.	Description	Packaging
Z31144168	Rack for I.V. pole, single - Devon™	1 ea/bx
Z5097	Rack with I.V. pole, double - Curity™	1 ea/ca
Z31144184	Rack for I.V. pole, triple - Devon™	1 ea/ca
Z31144176	Wall mounted rack, single - Devon™	1 ea/ca
Z31144192	Wall mounted rack, triple - Devon™	1 ea/ca

Transfer Products

Devon™ Blade Shield Scalpel Holder


- Sterile, disposable blade shield holds sharp edges of scalpels in the "down" position, helping to prevent accidental injuries
- Holds up to three scalpels in a convenient, ready-to-use position
- Adhesive tabs allow blade shields to be secured to any work surface
- Not made with natural rubber latex



Item No.	Description	Packaging
Z31142576	Blade Shield Scalpel Holder - Devon™	12 ea/bx, 8 bx/ca

Lap Sponges

OR Essentials





Z31142527

Anti-Fog Solution

- Anti-fog solution is sterile, non-toxic, laboratory tested and safe for use within body cavities
- May be used on endoscopic lenses, surgical cameras, protective eye wear, goggles or glasses
- Inhibits condensation and clouding
- Adhesive backed sponge and needle counter may be securely attached to any work surface
- Screw cap, individually wrapped, sterile
- Not made with natural rubber latex

Item No.	Description	Packaging
With foam pad		
Z31142527	6-gram anti-fog solution - Devon™	12 ea/bx, 4 bx/ca



Z31140299


Z31140315

Devon™ Surgical Clamps

- Disposable plastic, available sterile or non-sterile
- Tube occluding, disposable clamp will not cut fine tubes
- Individually wrapped
- Not made with natural rubber latex

Disposable Drape Clamps

Item No.	Description	Packaging
Z31140281	Drape clamp, plastic, sterile - Devon™	10 ea/bx, 8 bx/ca
Z31140299	Drape clamp, plastic, non-sterile, bulk - Devon™	100 ea/ca
Z31140315	Tube occluding clamp, plastic, non-sterile, bulk - Devon™	50 ea/bx, 1 bx/ca



Z31142790

Cautery Tip Cleaning Systems

- Cautery cleaning pad - safely and effectively cleans electrosurgical pencil tips, strong adhesive backing, radiopaque
- Micro tip wipe, lint and particulate-free foam wipe for cleaning the fine tips of micro instruments and lenses
- Sterile individually wrapped unless specified
- Not made with natural rubber latex

Item No.	Description	Packaging
300-255	Cautery cleaning pad	50 ea/bx, 4 bx/ca
Z31142790	Cautery cleaning pad - Devon™	25 ea/bx, 4 bx/ca
Z31142717	Cautery cleaning pad - Devon™	80 ea/bx, 4 bx/ca
Z31142535	Micro tip wipe, foam - Devon™	25 ea/bx, 4 bx/ca
Z31160099	Micro tip wipe, bulk, non-sterile	300 ea/ca



820-3409, 820-3609

Esmark Bandages

- Provides necessary compression for limb exsanguination during bloodless procedures
- Available in 4 and 6 in. widths
- Not made with natural rubber latex

Item No.	Description	Packaging
820-3409	4 in. x 3 yd., blue	20 ea/ca
820-3412	4 in. x 4 yd., blue	20 ea/ca
820-3609	6 in. x 3 yd., blue	20 ea/ca
820-3612	6 in. x 4 yd., blue	20 ea/ca



Insufflation Tubing

- Maintain pneumoperitoneum during abdominal endoscopy
- Filter capable of filtering particles down to 0.1 microns
- Flexible, lightweight tubing helps prevent kinking
- Sterile individually wrapped
- Not made with natural rubber latex

Item No.	Description	Packaging
C-2003-ST	10 ft. push-on fitting, with hi-flow filter	25 ea/ca
C-2005-ST	10 ft. quick connect fitting, with hi-flow filter	25 ea/ca
C-2025-ST	10 ft. quick connect fitting, with pressure sense line	25 ea/ca
C-2000-ST	10 ft. rotating luer lock, with hi-flow filter	25 ea/ca
C-2012-ST	12 ft. rotating luer lock, with hi-flow filter	25 ea/ca

Light Handle Systems

- Designed to function as an extension of customer specific OR lights. Adapter available for all popular light models
- Provides three choices to extend the sterile field: flexible, rigid, or semi-rigid light handle covers
- Convenient and cost-effective tool to drive standardization and minimize sterile reprocessing
- Not made with natural rubber latex



Universal Light Handles (ULH) and Universal Light Adapters (ULA)

- ULH is designed for use with flexible or semi-rigid light handle covers for Cardinal Health™ Brand Products only.
- ULA is designed for use with rigid light handles for Cardinal Health™ Brand Products only.

Adapters are customized to your light. Please contact your Cardinal Health sales representative for details

Flexible Light Handle Covers

- Sterile, soft cover slides over Universal Light Handle
- Single use only
- Not made with natural rubber latex



5128-FG

Item No.	Description	Packaging
5128-FG	Light handle cover, flexible (1 per pack)	1/pk, 32 pk/bx, 4 bx/ca
5160-2FG	Light handle cover, flexible (2 per pack)	2/pk, 20 pk/bx, 4 bx/ca

Rigid Light Handle

- Sterile, handle screws into Universal Light Adapter
- Single use only
- Not made with natural rubber latex



524-RG

Item No.	Description	Packaging
Z31140133	Rigid light handle (1 per pack) - Devon™	1/pk, 16 pk/bx, 8 bx/ca
Z31140125	Rigid light handle (2 per pack) - Devon™	2/pk, 14 pk/bx, 4 bx/ca
Z31140117	Rigid light handle (3 per pack) - Devon™	3/pk, 7 pk/bx, 4 bx/ca
CH524-RG	Rigid light handle (1 per pack)	1/pk, 24 pk/ca
CH548-2RG	Rigid light handle (2 per pack)	2/pk, 24 pk/ca

Semi-Rigid Light Sleeve

- Sterile, sturdy, flexible cover for Universal Light Handle
- Single use only
- Not made with natural rubber latex

Item No.	Description	Packaging
Z31140273	Semi-rigid light sleeve (1 per pack) - Devon™	1/pk, 70 pk/ca



250SM

Marking Pens, Surgical

Mini Marking Pens

- Intended for surgical marking on intact skin
- Non-toxic, non-smear, non-irritating Gentian violet ink
- Sterile and non-sterile versions available
- Not made with natural rubber latex

Item No.	Description	Packaging
Regular-tip		
250SM	Mini marking pen, sterile	50 ea/bx, 12 bx/ca
250SMNS	Mini marking pen, non-sterile	100 ea/bx, 12 bx/ca
Z31146010*	Mini marking pen, non-sterile - Devon™	100 ea/bx, 2 bx/ca

**Packaged in a cylindrical, clear plastic case for ease of pen selection.*

Prep Resistant Surgical Marking Pens

- ChloraPrep® resistant Gentian Violet ink increases visibility of surgical markings used for correct site surgery
- Ink specifically formulated for use with ChloraPrep® products, eliminating need for multiple applications
- Intended for surgical marking on intact skin
- Wider barrel than traditional marking pens
- Not made with natural rubber latex

Item No.	Description	Packaging
Mini-regular-tip		
Z31146020	Mini prep resistant surgical marking pen, non-sterile - Devon™	100 ea/bx, 2 bx/ca

See more marking pens next page



Z31145777

Marking Pens, Surgical continued

Traditional Surgical Marking Pens

- Non-toxic, non-smear, non-irritating Gentian violet ink
- Intended for surgical marking on intact skin
- Choice of regular, fine or tapered tip
 - Convex cone shape of tapered tip marking pens allows user to control and vary the width of markings
- Available with a flexible ruler, 9 blank labels or both ruler and labels
- Not made with natural rubber latex

Item No.	Description	Packaging
Dual-tip with ruler		
CH250P	Dual-tip marker, 4 cm barrel ruler	50 ea/bx, 12 bx/ca
Fine-tip		
Z31145942	Traditional marking pen, 9 cm ruler cap - Devon™	25 ea/bx, 4 bx/ca
Fine-tip with ruler and labels		
C250FPRL	Traditional marking pen, 6 in. barrel ruler	50 ea/bx, 12 bx/ca
Regular-tip		
Z31145819	Traditional marking pen - Devon™	25 ea/bx, 4 bx/ca
Z31160503	Traditional marking pen, non-sterile - Devon™	100 ea/ca
Z31145926	Traditional marking pen, 9 cm ruler cap - Devon™	25 ea/bx, 4 bx/ca
Regular-tip with ruler		
CH250GPR	Traditional marking pen	50 ea/bx, 12 bx/ca
Z31145785	Traditional marking pen - Devon™	25 ea/bx, 4 bx/ca
CH212WPRA	Traditional marking pen, 4 cm barrel ruler	12 ea/bx
Regular-tip with labels		
Z31145918	Traditional marking pen, 9 cm ruler cap - Devon™	25 ea/bx, 4 bx/ca
Z31160552	Traditional marking pen, 9 cm ruler cap, non-sterile - Devon™	100 ea/ca
Regular-tip with ruler and labels		
CH250GPRL	Traditional marking pen, 4 cm barrel ruler	50 ea/bx 12 bx/ca
Z31145777	Traditional marking pen, non-sterile - Devon™	25 ea/bx, 4 bx/ca
Taper-tip with ruler		
CH212PR	Traditional marking pen, 4 cm barrel ruler	12 ea/bx
CH250PR	Traditional marking pen, 4 cm barrel ruler	50 ea/bx, 12 bx/ca
Taper-tip with ruler and labels		
CH250PRL	Traditional marking pen, 4 cm barrel ruler	50 ea/bx, 12 bx/ca

See more marking pens next page



CH212MIL

Marking Pens, Surgical Dual Tip

Dual Tip: Surgical/Utility Marker

- Dual tip surgical marker include both regular and fine tips marking options in one pen.
- Cost effective method to decrease surgical waste
- Available with a flexible ruler, 9 blank labels or both ruler and labels
- Sterile, individually wrapped, unless otherwise noted
- Not made with natural rubber latex

Item No.	Description	Packaging
Z31145900	Dual tip marker, 9 cm ruler cap - Devon™	25 ea/bx, 4 bx/ca
With labels		
CH212MIL	Dual-tip marker, 4 cm barrel ruler	12 ea/bx, 56 bx/ca
Z31145868	Dual-tip marker, 9 cm ruler cap - Devon™	25 ea/bx, 4 bx/ca
With ruler		
Z31160560	Dual-tip marker, 9 cm ruler cap - Devon™	25 ea/bx, 4 bx/ca



Z31145884

Marking Pens, Utility

Utility Marking Pens

- May be used on most surfaces, including non-woven CSR wrap, sterilization tape, plastic, glass, or stainless steel
- Non-smearing, fast drying, non-toxic ink
- Non-toxic ink may be sterilized via steam or gas
- Available with a flexible ruler, 9 blank labels or both ruler and labels
- Sterile individually wrapped, unless otherwise noted
- Not made with natural rubber latex

Item No.	Description	Packaging
Fine-tip		
Z51000177	Utility marking pen	50 ea/ca
Fine-tip with labels		
Z31145884	Utility marking pen - Devon™	25 ea/bx, 4 bx/ca
Regular-tip with labels		
CH225UTL	Utility marking pen	25 ea/bx, 4 bx/ca
CHUT-50	Utility marking pen, non-sterile	50 ea/ca
Taper tip with labels		
CH225UTFL	Utility marking pen	25 ea/bx, 4 bx/ca



Medication Labels

- Medication labels enable healthcare providers to comply with patient safety goals
- We offer a portfolio of medication labels to meet a variety of provider needs.
- Universal labels for syringes, medication cups, bowls and basins
- Flexibility to choose specific labels for your department
- Waterproof so safety is not compromised

Item No.	Description	Packaging
ZCT6006-1	Chemotherapy caution label, 3 x 0.5 in.	500 ea/rl, 4 rl/ca
ZCT6011-1	Chemotherapy caution label, 3 x 1 in.	500 ea/rl, 4 rl/ca
SPP99IDGNA	Generic label sheet	50 ea/ca
SPP99MLGNA	Generic label sheet, with initial and date	50 ea/ca
Z31145769	Removable medication labels, 0.5 x 2 in., white, non-sterile	3000 ea/ca
SPP99MLCLA	Standard Cath Lab labeling sheet	50 ea/ca
SPP99IDCLA	Standard Cath Lab labeling sheet, with initials and date	50 ea/ca
SPP99MLCVA	Standard CV OR label sheet	50 ea/ca
SPP99IDENA	Standard ENT labeling sheet, with initials and date	50 ea/ca
SPP99MLEPA	Standard EP labeling sheet	50 ea/ca
SPP99MLOPA	Standard Ophthalmic label sheet	50 ea/ca
SPP99MLORA	Standard OR label sheet	50 ea/ca
SPP99IDORA	Standard OR label sheet, with initials and date	50 ea/ca
SPP99IDOTA	Standard Orthopedic label sheet, with initials and date	50 ea/ca



OR Plastics

- Choose from an array of shapes and capacities
- Not made with natural rubber latex
- Sterile and non-sterile

Item No.	Description	Packaging
Sterile		
SSK9005A	Emesis basin, 700 mL	30 ea/ca
SSK9012A	Medicine cup, 2 oz. with lids	100 ea/ca
SSK9004A	Pitcher, 1200 mL	30 ea/ca
Z61000	Solution bowl, 16 oz.	75 ea/ca
Z61200	Solution bowl, 32 oz.	50 ea/ca
SSK9001A	Sponge bowl, 16 oz.	33 ea/ca
SSK9002A	Sponge bowl, 32 oz.	30 ea/ca
Non-sterile		
Z31160206	Clear heart basin, 1500cc, non-sterile	225 ea/ca
Z65000	Solution bowl, 16 oz., non-sterile	500 ea/ca
Z65600	Solution bowl, 32 oz., non-sterile	250 ea/ca

OR Essentials



Universal Blade/Handle Sizing Guide

Blade size	Compatible handle size
11	3, 3L, 3LA, 5, 7, 9
15	3, 3L, 3LA, 5, 7, 9

Surgical Blades

- Sterile, single use surgical blade
- Individually wrapped; available loose or in boxes of 10 or 12
- Not made with natural rubber latex

Item No.	Description	Packaging
Disposable Safety Scalpel		
Z139091	#11 blade, bulk	100 ea/ca
Z139093	#15 blade, bulk	100 ea/ca

Surgical Vessel Loops

- Silicone, radiopaque
- Available in four colours and two lengths
- Not made with natural rubber latex
- Vessel loops packaged two per pouch loosely wound in a pouch
- Devon™ vessel loops packaged two per pouch, securely wound around a holding card

Mini Surgical Vessel Loops

Item No.	Description	Packaging
HT9945702A	Blue, 17 in. (2 per pack), bulk	2 ea/pk, 100 pk/ca
HT9945710A	Red, 17 in. (2 per pack), bulk	2 ea/pk, 100 pk/ca
HT9945660A	White, 17 in. (2 per pack), bulk	2 ea/pk, 100 pk/ca
HT9945728A	Yellow, 17 in. (2 per pack), bulk	2 ea/pk, 100 pk/ca
Z31145702	Blue, 18 in. (2 per pack) - Devon™	2 ea/pk, 10 pk/bx, 10 bx/ca
Z31145710	Red, 18 in. (2 per pack) - Devon™	2 ea/pk, 10 pk/bx, 10 bx/ca
Z31145728	White, 18 in. (2 per pack) - Devon™	2 ea/pk, 10 pk/bx, 10 bx/ca
Z31145660	Yellow, 18 in. (2 per pack) - Devon™	2 ea/pk, 10 pk/bx, 10 bx/ca

Maxi Surgical Vessel Loops

Item No.	Description	Packaging
HT9945678A	Blue, 17 in. (2 per pack), bulk	2 ea/pk, 100 pk/ca
HT9945744A	Red, 17 in. (2 per pack), bulk	2 ea/pk, 100 pk/ca
HT9945686A	White, 17 in. (2 per pack), bulk	2 ea/pk, 100 pk/ca
HT9945694A	Yellow, 17 in. (2 per pack), bulk	2 ea/pk, 100 pk/ca
Z31145678	Blue, 18 in. (2 per pack) - Devon™	2 ea/pk, 10 pk/bx, 10 bx/ca
Z31145744	Red, 18 in. (2 per pack) - Devon™	2 ea/pk, 10 pk/bx, 10 bx/ca
Z31145686	White, 18 in. (2 per pack) - Devon™	2 ea/pk, 10 pk/bx, 10 bx/ca
Z31145694	Yellow, 18 in. (2 per pack) - Devon™	2 ea/pk, 10 pk/bx, 10 bx/ca



Miscellaneous

Item No.	Description	Packaging
C28000-020	Rubber bands, sterile (2 per pack)	2/pk, 100 pk/cs

OR Lap Sponges and Towels

OR Lap Sponges and Towels



6545, 6556



4048

A28700-001 Blue, A28200-001 Green, 28300-002 White

Curity™ Lap Sponges

- 100% cotton disposable lap sponges
- X-Ray Detectable
- Tray or soft pack available
- Sterile and non-sterile

Item No.	Description	Packaging
Curity™ non-sterile, prewashed sponges in bulk pack		
Z6822	18 x 18 in., prewashed, non-sterile, loop	200 ea/ca
Curity™ sterile, prewashed sponges in soft pack		
Z6545	4 x 18 in., prewashed, sterile, loop	5 ea/pk, 40 pk/ca
Z6556	4 x 36 in., prewashed, sterile, loop	5 ea/pk, 40 pk/ca
Z6534	12 x 12 in., prewashed, sterile, loop	5 ea/pk, 40 pk/ca
Z6522	18 x 18 in., prewashed, sterile, loop	5 ea/pk, 40 pk/ca
Curity™ sterile, prewashed sponges in tray pack		
Z6045	4 x 18 in., prewashed, sterile, loop	5 ea/pk, 40 pk/ca
Z6056	8 x 36 in., prewashed, sterile, loop	5 ea/pk, 40 pk/ca
Z6034	12 x 12 in., prewashed, sterile, loop	5 ea/pk, 40 pk/ca
Z6022	18 x 18 in., prewashed, sterile, loop	5 ea/pk, 40 pk/ca

Kerlix™ Lap Sponges

- Dual-colored radiopaque handles that remain visible when saturated with blood
- Unique sizes
- Sterile only
- Tray pack only

Item No.	Description	Packaging
Kerlix™ Sterile, prewashed sponges in tray pack		
Z4048	3 x 17 in., prewashed, sterile, loops	5 ea/pk, 40 pk/ca
Z4037	8 x 32 in., prewashed, sterile, loops	5 ea/pk, 40 pk/ca
Z4031	17 x 26 in., prewashed, sterile, loops	5 ea/pk, 40 pk/ca

OR Towels

Standard OR Towels

OR Lap Sponges and Towels


- Available in blue, green and white
- Multiple packaging options available to help drive standardization

Item No.	Description	Packaging
Sterile blue		
A28700-001	OR towel 17 x 24 in., blue	1 ea/pk, 40 pk/ca
A28700-002	OR towel 17 x 24 in., blue	2 ea/pk, 40 pk/ca
A28700-004	OR towel 17 x 24 in., blue	4 ea/pk, 20 pk/ca
A28700-006	OR towel 17 x 24 in., blue	6 ea/pk, 12 pk/ca
Non-Sterile blue		
A28100-999	OR towel 17 x 28 in., blue	100 ea/ca
Sterile green		
A28200-001	OR towel 17 x 24 in., green	1 ea/pk, 40 pk/ca
A28200-002	OR towel 17 x 24 in., green	2 ea/pk, 40 pk/ca
A28200-004	OR towel 17 x 24 in., green	4 ea/pk, 20 pk/ca
Sterile white		
A7550	OR towel 15.5 x 23 in., white	2 ea/pk, 120 pk/ca
28300-002	OR towel 17 x 24 in., white	2 ea/pk, 40 pk/ca
A28300-004	OR towel 17 x 24 in., white	4 ea/pk, 20 pk/ca

X-Ray Detectable OR Towels

- Available in blue, green and white
- Radiopaque tag sewn into each towel

Item No.	Description	Packaging
Sterile blue		
28700-002X	X-Ray detectable towel, 17 x 27 in., blue	2 ea/pk, 40 pk/ca
28700-004X	X-Ray detectable towel, 17 x 24 in., blue	4 ea/pk, 20 pk/ca
28700-006X	X-Ray detectable towel, 17 x 24 in., blue	6 ea/pk, 12 pk/ca
Sterile green		
28200-004X	X-Ray detectable towel, 17 x 24 in., green	4 ea/pk, 20 pk/ca
Sterile white		
28300-002X	X-Ray detectable towel, 17 x 27 in., white	2 ea/pk, 40 pk/ca
28300-004X	X-Ray detectable towel, 17 x 24 in., white	4 ea/pk, 20 pk/ca
28300-006X	X-Ray detectable towel, 17 x 24 in., white	6 ea/pk, 12 pk/ca



SPP99UB2PK

Universal Low Lint Towels

- 99% less lint generated compared to cotton towels
- Absorbs up to 30% more fluid at equal or faster speed compared to cotton towels
- Handles and folds similar to cotton

Item No.	Description	Packaging
UB017026	Universal low lint towel, 17 x 26 in., blue, non-sterile	240 ea/ca
SPP99UB2PK	Universal low lint towel, 17 x 26 in., blue, sterile	2 ea/pk, 70 pk/ca
SPP99UB4PK	Universal low lint towel, 17 x 26 in., blue, sterile	4 ea/pk, 70 pk/ca

For more information or to order, contact your local
Cardinal Health Canada **sales representative today**
or **call 1.888.291.5033**