



# Chemotherapy Protection and Safety Products Catalogue



# Table of Contents

	<b>Protective Apparel</b> ..... 4
	ChemoPlus™ Chemotherapy Gowns ChemoPlus™ Chemotherapy Sleeves Cardinal Health™ Chemotherapy Gowns ChemoPlus™ Protective Eyewear
	<b>Nitrile Exam Gloves</b> ..... 6
	Cardinal Health™ Decontamination Nitrile Exam Gloves Cardinal Health™ Extended Cuff Nitrile Exam Gloves
	<b>Absorbent Products</b> ..... 7
	ChemoPlus™ Absorbent Pads ChemoPlus™ Low Lint Towels ChemoPlus™ Prep Mats ChemoPlus™ Floor Mats
	<b>Waste Containers, Bags and Labels</b> ..... 8
	ChemoPlus™ Chemo Soft Waste Container ChemoPlus™ Waste Bags ChemoPlus™ Drug Transport Bags ChemoPlus™ Medical Labels
	<b>Chemotherapy Sharps Containers</b> ..... 10
	ChemoSafety™ Sharps Container with Slide Lid or Hinged Lid
	<b>IVA Security Seals</b> ..... 11
	ChemoPlus™ IVA Seals for Bottles & Vials ChemoPlus™ IVA Seals for I.V. Bags ChemoPlus™ IVA Seals for Syringes & Medication Containers
	<b>Vial Venting System</b> ..... 13
	ChemoPlus™ Vial Vent Chemo Dispensing Pin
	<b>Index</b> ..... 14

CHEMOTHERAPY PROTECTION THAT WORKS AS HARD AS YOU DO. Healthcare professionals work hard to **protect** and **treat patients**. They deserve chemotherapy protection that works hard to **protect them**.

As cancer treatment evolves, so do the drug touch points. Exposure risks exist in many areas of a healthcare facility.

Our comprehensive line of chemotherapy management products offers healthcare professionals with protection from potential hazards associated with exposure to chemotherapy drugs at every step of care; from preparation and handling to administration and disposal.



# Protect yourself and others



## Know your exposure risk

# 75,000

Canadians are estimated as occupationally exposed to antineoplastic agents<sup>1</sup>, putting more healthcare workers at risk to hazardous drug exposure<sup>1</sup>.

### Short-term risks<sup>2</sup>

- Allergic asthma
- Ocular irritation
- Headache
- Cough
- Dizziness
- Nausea/vomiting
- Hair thinning
- Skin rash

### Long-term risks<sup>2</sup>

- DNA damage
- Increased cancer risk
- Infertility
- Miscarriage
- Birth defects
- Organ damage

## How and where contamination occurs



### Hazardous drugs

**More than half** of hazardous drugs are antineoplastic agents used in cancer treatment, which are among the **most toxic** administered in healthcare settings<sup>3</sup>.



### Dermal exposure

For healthcare workers, dermal exposure —whether direct or indirect — is the **most common** route of contaminant entry<sup>3</sup>.



### Bodily fluids

Patients can excrete chemo drugs for up to **7 days** after administration<sup>3</sup>; healthcare workers caring for them are **repeatedly exposed** to bodily fluids<sup>4</sup>.



### High-touch surfaces

Hazardous drugs are found on up to **61%** of surfaces<sup>3</sup>, including:

- **Shared surfaces** touched by healthcare workers, patients and family members
- **Patient room surfaces**, including toilet seats, doorknobs and floors



**Smart protection at every step of care**

**Potential exposure to hazardous drugs can occur in every aspect if proper precautions are not taken<sup>5,6</sup>:**

- Receipt, storage and transport
- Drug compounding
- Administration
- Handling body fluids post-administration
- Deactivation and/or decontamination
- Cleaning and disinfecting
- Spill management
- Waste disposal

To order, contact your sales representative or visit [CardinalHealth.ca/ChemoPlusGowns](https://www.cardinalhealth.ca/ChemoPlusGowns)

1. Estimating National-Level Exposure to Antineoplastic Agents in the Workplace: CAREX Canada Findings and Future Research Needs. <https://academic.oup.com/annweh/article/61/6/656/3863057> Submitted 20 January, 2017; revised 29 April, 2017; editorial decision 1 May 2017; revised version accepted 11 May 2017.

2. Olsen MM, et al. *Safe Handling of Hazardous Drugs*. Pittsburgh, PA: ONS; 2018.

3. Walton AL, et al. *Oncol Nurs Forum*. 2020;47(3):263-272.

4. Walton AL, et al. *Oncol Nurs Forum*. 2019;46(6):689-700.

5. Antineoplastic agents: hazardous drug exposures in healthcare. CDC. <https://www.cdc.gov/niosh/topics/hazdrug/antineoplastic.html>. Published September 13, 2017. Accessed December 6, 2021

6. USP General Chapter <800> Hazardous Drugs - Handling in Healthcare Settings <https://www.usp.org/usp-chapter-800-download>

## Protective Apparel

• Cardinal Health ChemoPlus™ Chemotherapy gowns are engineered to stay at the forefront of regulatory standards and industry best practices

### ChemoPlus™ Chemotherapy Gowns

- ChemoPlus™ poly-coated chemotherapy apparel in light blue
- Full back style gowns that close and tie in the back and sleeve options
- Open back gowns that provide front coverage and breathability
- AAMI Level 3 gowns and sleeves meet strict standards for moderate fluid protection
- One of the first infection control apparel to offer sterile packaging options
- Evaluated for resistance to permeation of 16 chemotherapy drugs per ASTM F739-20<sup>1</sup>
- Meets USP <800> standards for handling hazardous drugs
- Found to be non-cytotoxic, non-irritating and non-sensitizing after being tested for biological clearance per ISO 10993

Item Code	Description	Size	Colour	Packaging
<b>Non-sterile, Full Back</b>				
ZCT5502T	ChemoPlus™ Chemotherapy Gowns	Medium	Light blue	11 ea/bg, 2 bg/ca
ZCT5503T	ChemoPlus™ Chemotherapy Gowns	Large	Light blue	11 ea/bg, 2 bg/ca
ZCT5504T	ChemoPlus™ Chemotherapy Gowns	X-Large	Light blue	11 ea/bg, 2 bg/ca
ZCT5505T	ChemoPlus™ Chemotherapy Gowns	XX-Large	Light blue	11 ea/bg, 2 bg/ca
<b>Non-sterile, Open Back</b>				
ZDP5003GT	ChemoPlus™ Chemotherapy Gowns	Medium	Light blue	14 ea/bg, 2 bg/ca
ZDP5001GT	ChemoPlus™ Chemotherapy Gowns	Large	Light blue	14 ea/bg, 2 bg/ca
ZDP5002GT	ChemoPlus™ Chemotherapy Gowns	X-Large	Light blue	14 ea/bg, 2 bg/ca
ZDP5004GT	ChemoPlus™ Chemotherapy Gowns	XX-Large	Light blue	14 ea/bg, 2 bg/ca
<b>Sterile, Full Back</b>				
ZCT5502TS	ChemoPlus™ Chemotherapy Gowns	Medium	Light blue	1 ea/bg, 14 bg/ca
ZCT5503TS	ChemoPlus™ Chemotherapy Gowns	Large	Light blue	1 ea/bg, 14 bg/ca
ZCT5504TS	ChemoPlus™ Chemotherapy Gowns	X-Large	Light blue	1 ea/bg, 14 bg/ca
ZCT5505TS	ChemoPlus™ Chemotherapy Gowns	XX-Large	Light blue	1 ea/bg, 14 bg/ca



ZCT5502T



ZDP5003GT

### ChemoPlus™ Universal Sleeve

- Evaluated for resistance to permeation of 16 chemotherapy drugs per ASTM F739-20<sup>1</sup>
- AAMI Level 3 sleeves meet strict standards for moderate fluid protection
- Found to be non-cytotoxic, non-irritating and non-sensitizing after being tested for biological clearance per ISO 10993

Item Code	Description	Size	Colour	Packaging
ZCT5500T	ChemoPlus™ Universal Sleeve, Non-Sterile	Universal	Blue	120 sleeves/ca
ZCT5500TS	ChemoPlus™ Universal Sleeve, Sterile	Universal	Blue	2 ea/bg, 60 bg/ca



ZCT5500T

 For further information please view our: [Chemotherapy Gowns Brochure](#)

## Poly-coated SMS Chemotherapy Gowns

- Engineered to provide users with protection from the potential hazards associated with handling chemotherapy agents
- Poly-coated material for durability and protection
- Lightweight, closed-back design, secured with ties for maximum comfort
- Low-linting material helps to provide a clean environment
- Long sleeves with knit cuffs provide added comfort
- Evaluated for resistance to permeation of 11 chemotherapy drugs per ASTM F739-12<sup>2</sup>



Item Code	Description	Size	Colour	Packaging
A8200CG	Poly-Coated SMS Chemotherapy Gowns	Universal	Blue	10 ea/pk, 10 pk/ca
A8201CG	Poly-Coated SMS Chemotherapy Gowns	X-Large	Blue	10 ea/pk, 10 pk/ca

## Cardinal Health™ Molded N95 Respirator/Surgical Mask

- NIOSH and FDA cleared for use in the OR during surgical procedures
- Meets CDC (Center for Disease Control and Prevention) guidelines for protection against TB
- Malleable aluminum nose wire
- Nylon/spandex elastic securement
- Not made with natural rubber latex
- Color White



N95A-ML

Item Code	Description	Size	Colour	Packaging
N95A-ML	Mask N95 STANDARD	Med-Large	White	20 ea/bx, 6 bx/ca

## ChemoPlus™ Protective Eyewear

- Splash resistant
- Universal fit
- Comfortably fits over most eyeglasses
- Ventilated



ZDP5030G

Item Code	Description	Size	Packaging
ZDP5030G	Protective Wrap-Around Goggles	Universal	6 ea/ca
ZCT0400-1	Plastic Eyewear with Side Shields	Universal	10 ea/bx, 5 bx/ca



ZCT0400-1

 For further information please view our: [Chemotherapy Gowns Brochure](#)

1. ASTM International. Standard Test Method for Permeation of Liquids and Gases Through Protective Clothing Materials Under Conditions of Continuous Contact. <https://www.astm.org/Standards/F739.htm>. Accessed October 6, 2021

2. ASTM International. Standard Test Method for Permeation of Liquids and Gases Through Protective Clothing Materials Under Conditions of Continuous Contact. <https://www.astm.org/Standards/F739.htm>. Accessed January 14, 2014

## Nitrile Exam Gloves

- Tested per ASTM D6978 for Use with Certain Chemotherapy Drugs and Tested for use with various lab chemicals<sup>1</sup>
- Compliant with USP <800> standards<sup>1</sup>

### Cardinal Health™ Decontamination Nitrile Exam Gloves

- Nominal Fingertip Thickness: 11 mil
- Full Finger Textured
- 16 in. cuff for added coverage and protection
- Powder-free
- Ambidextrous
- Not made with natural rubber latex
- Tested for Puncture Resistance<sup>2</sup>



88NDS

Item Code	Description	Size	Thickness	Packaging
88NDS	Decontamination Nitrile Exam Gloves, Blue	Small	11 mil	50 ea/bx, 10 bx/ca
88NDM	Decontamination Nitrile Exam Gloves, Blue	Medium	11 mil	50 ea/bx, 10 bx/ca
88NDL	Decontamination Nitrile Exam Gloves, Blue	Large	11 mil	50 ea/bx, 10 bx/ca
88NDXL	Decontamination Nitrile Exam Gloves, Blue	X-Large	11 mil	50 ea/bx, 10 bx/ca

### Cardinal Health™ Extended Cuff Nitrile Exam Gloves

- Nominal Fingertip Thickness: 6.9 mil
- Fingertip Textured
- 12 in. cuff for added coverage and protection
- Ambidextrous
- Not made with natural rubber latex



N8851XPB

Item Code	Description	Size	Thickness	Packaging
N8851XPB	Extended Cuff Nitrile Exam Gloves, Blue	Small	6.9 mil	100 ea/bx, 10 bx/ca
N8852XPB	Extended Cuff Nitrile Exam Gloves, Blue	Medium	6.9 mil	100 ea/bx, 10 bx/ca
N8853XPB	Extended Cuff Nitrile Exam Gloves, Blue	Large	6.9 mil	100 ea/bx, 10 bx/ca
N8854XPB	Extended Cuff Nitrile Exam Gloves, Blue	X-Large	6.9 mil	100 ea/bx, 10 bx/ca

<sup>1</sup> Meets the ASTM D6978 testing standard for resistance to permeation of various chemotherapy drugs. This Cardinal Health product helps healthcare facilities comply with the USP <800> standard for the safe handling of hazardous drugs.

<sup>2</sup> This glove has been tested for puncture resistance per ASTM F1342, Standard Test Method for Protective Clothing Material Resistance to Puncture.

## Absorbent Products

### ChemoPlus™ Absorbent Pads

- Highly absorbent
- Pads are impregnated with Green-Z™ solidifying agent to instantly turn liquid spills into a semi-gel
- Latex free, single use




ZBS0013

Item Code	Description	Size	Packaging
ZBS0013	Absorbent Pad	12 x 12 in. (30.5 x 30.5 cm)	100 ea/ca

### ChemoPlus™ Low Lint Towels

Item Code	Description	Size	Packaging
ZCT0014	Low Lint Towel	9 x 9 in. (23 x 23 cm)	300 ea/ca



ZCT0014

### ChemoPlus™ Prep Mats

- Highly absorbent design
- Available in several sizes
- Economy version available
- Latex free




ZCT0300-1

Item Code	Description	Size	Packaging
ZDP5010M	Economy Absorbent Prep Mat	17 x 19 in. (43.2 x 48.3 cm)	125 ea/ca
ZCT0300-1	Prep Mat Spill Absorbent	11 x 17 in. (27.9 x 43.2 cm)	50 ea/ca, 5 bx/ca
ZCT0302-1	Prep Mat Spill Absorbent	16 x 22 in. (40.6 x 55.9 cm)	50 ea/ca, 5 bx/ca
ZCT0303-1	Prep Mat Spill Absorbent	16 x 69 in. (40.6 cm x 1.8 m)	15 ea/ca, 5 bx/ca

### ChemoPlus™ Floor Mats

- Convenient and economical way to trap impurities from foot traffic in the pharmacy
- Easy to adhere to floor
- Simple to peel back, exposing new layer
- Thirty peel-back layers per mat
- Latex free



ZCT0007

Item Code	Description	Size	Colour	Packaging
ZCT0007	Adhesive Floor Mat	18 x 46 in. (45.7 cm x 1.2 m)	White	4 ea/ca
ZCT0071	Adhesive Floor Mat	18 x 46 in. (45.7 cm x 1.2 m)	Blue	4 ea/ca

## Waste Containers, Bags and Labels

### ChemoPlus™ Chemo Soft Waste Container

- Easy to assemble
- Dual flap opening with self-locking lid
- Holds approximately 20 gal. of waste
- ZCT1500 includes one 20 gal. waste bag
- Latex free

Item Code	Description	Size	Thickness	Colour	Packaging
ZCT1500	Chemo Soft Waste Container	20 gal.	4 mil	Yellow	6 ea/ca



ZCT1500

### ChemoPlus™ Waste Bags

- Available in a variety of sizes
- 4 mil thick, very durable
- Chemotherapy yellow color
- Closure ties included
- Latex free

Item Code	Description	Size	Thickness	Colour	Packaging
ZDP5041B	Chemo Soft Waste Bag	2 gal, 14 x 20 in.	4 mil	Yellow	250 ea/ca
ZDP5042B	Chemo Soft Waste Bag	15 gal, 18 x 6 x 33 in.	4 mil	Yellow	100 ea/ca
ZCT2100	Chemo Soft Waste Bag	20 gal, 23 x 8 x 41 in.	4 mil	Yellow	100 ea/ca
ZDP5043B	Chemo Soft Waste Bag	30 gal, 20 x 13 x 44 in.	4 mil	Yellow	100 ea/ca
WCHE31x43	Chemo Soft Waste Bag	31 x 41 in.	3 mil	Yellow	100 ea/ca



ZCT2100

### ChemoPlus™ Drug Transport Bags

- Available in two sizes
- 4 mil thick for maximum protection
- Zip-lock closure
- Clear to visualize contents
- Latex free

Item Code	Description	Size	Thickness	Colour	Packaging
ZCT0575	Chemo Drug Transport Bag	6 x 9 in. (15.2 x 22.9 cm)	4 mil	White	200 ea/ca
ZCT0500	Chemo Drug Transport Bag	12 x 15 in. (30.5 x 38.1 cm)	4 mil	White	250 ea/ca



ZCT0575



## ChemoPlus™ Medical Labels

- Includes labels for pharmacy, oncology, I.V. therapy, and patient care
- Reduces risk of medication errors and ensures clinician and patient safety
- Pre-printed message provides convenient and cost-effective communication
- "Caution: Chemotherapy Handle with gloves, Dispose of properly"
- Latex free



ZCT6011-1

Item Code	Description	Size	Packaging
ZCT6006-1	Chemo Caution Label	1/2 x 3 in. (1.3 x 7.6 cm)	500 ea/bx, 4 bx/ca
ZCT6011-1	Chemo Caution Label	1 x 3 in. (2.5 x 7.6 cm)	500 ea/bx, 4 bx/ca

# Chemotherapy Sharps Containers

## ChemoSafety™ Sharps Containers with Slide Lid or Hinged Lid

- For the disposal of trace chemotherapy waste
- Sealing gasket lids provide leak-resistant containment of fluid
- Absorbent pad in bottom of container absorbs excess liquid
- Adjustable sliding lid reduces the possibility of aerosolization and accommodates chemotherapy waste of various types and sizes
- Available carts with wheels and built-in carrying handles for easy transport and lifting of heavy, full containers
- Containers provide temporary and final locking closure for disposal
- Nestable containers save valuable storage space
- Sliding gasket lids have a round 6 in. opening to accommodate up to 1000 mL bottles
- Bilingual label (French, English)



Item Code	Description	Size (H x D x W)	Colour	Packaging
<b>Sliding Lid</b>				
Z8985CRS	Sharps Container	8 gal. (17.75 x 11 x 15.5 in.)	Red	10 ea/ca
Z8934CRS	Sharps Container	12 gal. (18.75 x 12.75 x 18.5 in.)	Red	10 ea/ca
Z8939CRS	Sharps Container	18 gal. (26 x 12.75 x 18.25 in.)	Red	5 ea/ca
<b>Hinge Lid</b>				
Z8982CRH	Sharps Container	2 gal. (10 x 7.25 x 10.5 in.)	Red	20 ea/ca
Z8985CRH	Sharps Container	8 gal. (17.75 x 11 x 15.5 in.)	Red	10 ea/ca
Z8989CRH	Sharps Container	18 gal. (26 x 12.75 x 18.25 in.)	Red	5 ea/ca

# IVA Security Seals

## ChemoPlus™ IVA Seals for Bottles & Vials

- Dual warning system: Seals self-destruct upon attempted removal and needle punctures are highly visible in laminated foil
- Sterile IVA seals protect point of entry of IV admixtures from contamination and accidental double dosing
- Multiple sizes and colours available to fit most bottles and vials
- Latex free & PVC free



Item Code	Description	Packaging
<b>For Top Vials</b>		
ZCP3005A	Silver IVA Seals, 9 mm	1100 ea/bx, 12 bx/ca
ZCP3005R	Red IVA Seals, 9 mm	1100 ea/bx, 12 bx/ca
ZCP3003A	Silver IVA Seals, 13 mm	1100 ea/bx, 20 bx/ca
ZCP3003B	Blue IVA Seals, 13 mm	1100 ea/bx, 20 bx/ca
ZCP3003G	Green IVA Seals, 13 mm	1100 ea/bx, 20 bx/ca
ZCP3003R	Red IVA Seals, 13 mm	1100 ea/bx, 20 bx/ca
ZCP3003Y	Yellow IVA Seals, 13 mm	1100 ea/bx, 20 bx/ca
ZCP3002G	Green IVA Seals, 20 mm	1000 ea/bx, 12 bx/ca
ZCP3002B	Blue IVA Seals, 20 mm	1000 ea/bx, 12 bx/ca
ZCP3002Y	Yellow IVA Seals, 20 mm	1000 ea/bx, 12 bx/ca
ZCP3002O	Orange IVA Seals, 20 mm	1000 ea/bx, 12 bx/ca
ZCP3002R	Red IVA Seals, 20 mm	1000 ea/bx, 12 bx/ca
<b>For Top Bottles &amp; Piggyback Containers</b>		
ZCP3001A	Silver IVA Seals, 28 mm	1000 ea/bx, 12 bx/ca
ZCP3001G	Green IVA Seals, 28 mm	1000 ea/bx, 12 bx/ca
ZCP3001B	Blue IVA Seals, 28 mm	1000 ea/bx, 12 bx/ca
ZCP3001R	Red IVA Seals, 28 mm	1000 ea/bx, 12 bx/ca
ZCP3001Y	Yellow IVA Seals, 28mm	1000 ea/bx 12 bx/ca
<b>For Top Bottles &amp; Lily's Faspak Containers</b>		
ZCP3009A	Silver IVA Seals, 36 mm	1000 ea/bx, 12 bx/ca
ZCP3009B	Blue IVA Seals, 36mm	1000 ea/bx, 12bx/ca
ZCP3009G	Green IVA Seals, 36 mm	1000 ea/bx, 12 bx/ca
ZCP3009R	Red IVA Seals, 36 mm	1000 ea/bx, 12 bx/ca

Most IVA seals are available in various colors:



A-Silver



R-Red



B-Blue



G-Green



O-Orange



Y-Yellow

## ChemoPlus™ IVA Seals for I.V. Bags

- Visual communication that prevents dosing errors
- Dual warning system: Seals self-destruct upon attempted removal and needle punctures are highly visible in laminated foil
- IVA seals on bags instantly alert staff that a drug has been added
- Sterile IVA seals protect point of entry of IV admixtures from contamination and accidental double dosing
- Three styles to accommodate most commonly used IV bag
- PVC free
- Latex free



Item Code	Description	Packaging
<b>For Hospira's Large Bags &amp; Add-Vantage™ System</b>		
ZCP3004A	Silver IVA Seals	1000 ea/bx, 12 bx/ca
ZCP3004B	Blue IVA Seals	1000 ea/bx, 12 bx/ca
ZCP3004G	Green IVA Seals	1000 ea/bx, 12 bx/ca
ZCP3004R	Red IVA Seal	1000 ea/bx, 12 bx/ca
<b>For Baxter's Viaflex™ &amp; Mini-Bag Plus Containers</b>		
ZCP3011A	Silver IVA Seals	1000 ea/bx, 20 bx/ca
ZCP3011G	Green IVA Seals	1000 ea/bx, 20 bx/ca
ZCP3011B	Blue IVA Seals	1000 ea/bx, 20 bx/ca
ZCP3011R	Red IVA Seals	1000 ea/bx, 20 bx/ca
ZCP3011Y	Yellow IVA Seals	1000 ea/bx, 20 bx/ca
<b>For B. Braun's Excel™ Bags &amp; Abbott's Small Bags</b>		
ZCP3013A	Silver IVA Seals	1000 ea/bx, 12 bx/ca
ZCP3013G	Green IVA Seals	1000 ea/bx, 12 bx/ca
ZCP3013B	Blue IVA Seals	1000 ea/bx, 12 bx/ca
ZCP3013R	Red IVA Seals	1000 ea/bx, 12 bx/ca
<b>For Baxter's Viaflex™ Bags &amp; Mini-Bag Plus Containers</b>		
ZCP3111Y	"Chemotherapy" Yellow/Black Ink	1000 ea/bx, 20 bx/ca



ZCT3004B



ZCT3013R

## ChemoPlus™ IVA Seals for Syringes & Medication Containers

- Ensures the security of needle covers, syringe caps, box tops and many other medication containers
- Non-sterile, tamper evident seals self-destruct upon attempted removal
- Unused medications may be returned and reissued with confidence if applied seals are intact
- Available in large, 4-1/4 in. (10.8 cm) and small, 3-5/8 in. (9.2 cm) lengths



Item Code	Description	Packaging
<b>For Larger Syringes, 4-1/4 in. length</b>		
ZCP3020B	Blue IVA Seals	1000 ea/bx, 12 bx/ca
ZCP3020Y	Yellow IVA Seals	1000 ea/bx, 12 bx/ca
ZCP3020R	Red IVA Seals	1000 ea/bx, 12 bx/ca
<b>For Smaller Syringes, 3-5/8 in. length</b>		
ZCP3030R	Red IVA Seals	1000 ea/bx, 12 bx/ca
<b>IVA™ for Smaller Syringes Seals</b>		
ZCP3730A	"Don't Use if Seal Is Broken" Silver/Red Ink	12 rl/ca (1000 ea/rl)
<b>IVA™ for Large Syringes Seals</b>		
ZCP3720A	"Don't Use if Seal Is Broken" Silver/Red Ink	1000 ea/bx, 12 bx/ca



ZCP3020B

## Vial Venting System

A system designed to prevent aerosolization of chemotherapy and other hazardous drugs into the atmosphere during reconstitution.

### ChemoPlus™ Vial Vent Chemo Dispensing Pin

- Can be used for reconstituting chemotherapy drugs, antibiotics, and other hazardous drugs
- 0.22 micron hydrophobic filter to safely remove any aerosolized drugs
- Safe, needleless system prevents needle stick injury
- Unique cannula is easy to insert
- Facilitates rapid transfer of solution
- Sterile, non-pyrogenic fluid path
- Luer-lock connection assures secure, safe connection to syringe
- Single use only
- Latex free



Item Code	Description	Packaging
ZDP5201V	Vial Vent Chemo Dispensing Pin	50 ea/bx, 4 bx/ca

# Index

Item Code	Product	Page	Item Code	Product	Page
88NDL	Decontamination Nitrile Exam Gloves, Large	6	ZCP3009A	Silver IVA Seals, 36 mm	11
88NDM	Decontamination Nitrile Exam Gloves, Medium	6	ZCP3009B	Blue IVA Seals, 36mm	11
88NDS	Decontamination Nitrile Exam Gloves, Small	6	ZCP3009G	Green IVA Seals, 36 mm	11
88NDXL	Decontamination Nitrile Exam Gloves, X Large	6	ZCP3009R	Red IVA Seals, 36 mm	11
A8200CG	Poly-Coated SMS Chemotherapy Gowns, Universal	5	ZCP3011A	Silver IVA Seals	12
A8201CG	Poly-Coated SMS Chemotherapy Gowns, Universal	5	ZCP3011B	Blue IVA Seals	12
N8851XPB	Extended Cuff Nitrile Exam Gloves. Small	6	ZCP3011G	Green IVA Seals	12
N8852XPB	Extended Cuff Nitrile Exam Gloves. Medium	6	ZCP3011R	Red IVA Seals	12
N8853XPB	Extended Cuff Nitrile Exam Gloves. Large	6	ZCP3011Y	Yellow IVA Seals	12
N8854XPB	Extended Cuff Nitrile Exam Gloves. X Large	6	ZCP3013A	Silver IVA Seals	12
N95A-ML	MASK N95 STANDARD	5	ZCP3013B	Blue IVA Seals	12
Z8934CRS	Sharps Container, Sliding Lid, 12 gal.	10	ZCP3013G	Green IVA Seals	12
Z8939CRS	Sharps Container, Sliding Lid, 18 gal.	10	ZCP3013R	Red IVA Seals	12
Z8982CRH	Sharps Container, Hinge Lid, 2 gal.	10	ZCP3020B	Blue IVA Seals	13
Z8985CRH	Sharps Container, Hinge Lid, 8 gal.	10	ZCP3020R	Red IVA Seals	13
Z8985CRS	Sharps Container, Sliding Lid, 8 gal.	10	ZCP3020Y	Yellow IVA Seals	13
Z8989CRH	Sharps Container, Hinge Lid, 18 gal.	10	ZCP3030R	Red IVA Seals	13
ZBS0013	Absorbent Pad	7	ZCP3111Y	Chemotherapy" Yellow/Black Ink	12
ZCP3001A	Silver IVA Seals, 28 mm	11	ZCP3720A	Don't Use if Seal Is Broken" Silver/Red Ink	13
ZCP3001B	Blue IVA Seals, 28 mm	11	ZCP3730A	Don't Use if Seal Is Broken" Silver/Red Ink	13
ZCP3001G	Green IVA Seals, 28 mm	11	ZCT0007	Adhesive Floor Mat, White	7
ZCP3001R	Red IVA Seals, 28 mm	11	ZCT0014	Low Lint Towel	7
ZCP3001Y	Yellow IVA Seals, 28mm	11	ZCT0071	Adhesive Floor Mat, Blue	7
ZCP3002B	Blue IVA Seals, 20 mm	11	ZCT0300-1	Prep Mat Spill Absorbent, 11 x 17 in.	7
ZCP3002G	Green IVA Seals, 20 mm	11	ZCT0302-1	Prep Mat Spill Absorbent, 16 x 22 in.	7
ZCP3002O	Orange IVA Seals, 20 mm	11	ZCT0303-1	Prep Mat Spill Absorbent, 16 x 69 in.	7
ZCP3002R	Red IVA Seals, 20 mm	11	ZCT0400-1	Plastic Eyewear with Side Shields	5
ZCP3002Y	Yellow IVA Seals, 20 mm	11	ZCT0500	Chemo Drug Transport Bag, 12 x 15 in.	8
ZCP3003A	Silver IVA Seals, 13 mm	11	ZCT0575	Chemo Drug Transport Bag, 6 x 9 in.	8
ZCP3003B	Blue IVA Seals, 13 mm	11	ZCT1500	Chemo Soft Waste Container, 20 gal.	8
ZCP3003G	Green IVA Seals, 13 mm	11	ZCT2100	Chemo Soft Waste Bag, 20 gal.	8
ZCP3003R	Red IVA Seals, 13 mm	11	ZCT5500T	ChemoPlus™ Universal Sleeve, Non-Sterile, Universal	4
ZCP3003Y	Yellow IVA Seals, 13 mm	11	ZCT5500TS	ChemoPlus™ Universal Sleeve, Sterile, Universal	4
ZCP3004A	Silver IVA Seals	12	ZCT5502T	ChemoPlus™ Chemotherapy Gowns, Medium	4
ZCP3004B	Blue IVA Seals	12	ZCT5502TS	ChemoPlus™ Chemotherapy Gowns, Medium	4
ZCP3004G	Green IVA Seals	12	ZCT5503T	ChemoPlus™ Chemotherapy Gowns, Large	4
ZCP3004R	Red IVA Seals	12	ZCT5503TS	ChemoPlus™ Chemotherapy Gowns, Large	4
ZCP3005A	Silver IVA Seals, 9 mm	11	ZCT5504T	ChemoPlus™ Chemotherapy Gowns, X-Large	4
ZCP3005R	Red IVA Seals, 9 mm	11	ZCT5504TS	ChemoPlus™ Chemotherapy Gowns, X-Large	4

# Index

<b>Item Code</b>	<b>Product</b>	<b>Page</b>
ZCT5505T	ChemoPlus™ Chemotherapy Gowns, XX-Large	4
ZCT5505TS	ChemoPlus™ Chemotherapy Gowns, XX-Large	4
ZCT6006-1	Chemo Caution Label, 1/2 x 3 in.	9
ZCT6011-1	Chemo Caution Label, 1 x 3 in.	9
ZDP5001GT	ChemoPlus™ Chemotherapy Gowns, Large	4
ZDP5002GT	ChemoPlus™ Chemotherapy Gowns, X-Large	4
ZDP5003GT	ChemoPlus™ Chemotherapy Gowns, Medium	4
ZDP5004GT	ChemoPlus™ Chemotherapy Gowns, XX-Large	4
ZDP5010M	Economy Absorbent Prep Mat	8
ZDP5030G	Protective Wrap-Around Goggles	5
ZDP5041B	Chemo Soft Waste Bag, 2 gal.	8
ZDP5042B	Chemo Soft Waste Bag, 15 gal.	8
ZDP5043B	Chemo Soft Waste Bag, 30 gal.	8
ZDP5201V	Vial Vent Chemo Dispensing Pin	13

For more information or to order, contact your local  
Cardinal Health Canada **sales representative today or call:**  
**Acute** - 1.888.291.5033  
**Continuing Care** - 1.877.878.7778

