

Salem Sump™ with ENFit™ Connection



The new **Salem Sump™** with an attached **ENFit™ connection** now has an extended indication for **suction and feeding**. Hospitals that have adapted to ENFit™ can now fully convert their practice. The new feeding indication makes this the **first ENFit™ compatible dual lumen suction/feeding tube in Canada**.

Product Features

- Suitable for gastric decompression, feeding, hydration and medication
- Gastric reflux prevention valve
- Radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Funnel end of vent lumen can be used to cap 5-in-1 adapter
- Centimeter markings
- Sterile



ENFit™ Connector

Enables use of the Salem Sump with ENFit™ connection system to keep patients safe.



Anti-Reflux Valve*

The anti-reflux valve helps prevent stomach content from exiting the vent lumen to keep clinicians safe.



Dual Lumen Tubing

The double lumen tube compensates for the negative pressure created by suction. It can now also be used for feeding.

Ordering Information

Salem Sump™ with ENFit™ Connection

Item No.	Description	Packaging
Z8888268060E	6 fr, length: 61 cm	10 ea/ca
Z8888268086E	8 fr, length: 61 cm	10 ea/ca
Z8888264911E	10 fr, length: 91 cm	50 ea/ca
Z8888264929E	12 fr, length: 122 cm	50 ea/ca

Item No.	Description	Packaging
Z8888264945E	14 fr, length: 122 cm	50 ea/ca
Z8888264960E	16 fr, length: 122 cm	50 ea/ca
Z8888264986E	18 fr, length: 122 cm	50 ea/ca
Z8888266197	Anti-Reflux Valve	10 ea/ca

*Anti-Reflux Valve is sold separately.

For more information or to order, please contact your Cardinal Health Canada sales representative or call 1.888.291.5033.