Wound Care. Anywhere.™

Cardinal Health™ NPWT Disposables and accessories

Our dressing kits and accessories have you covered regardless of which Cardinal Health™ NPWT system you choose.

Wound dressing kits

Black foam kit — open cell, reticulated, hydrophobic, high tensile strength

- Small, medium, large, and extra-large to handle a wide range of wound types and sizes
- Polyurethane Drape — easy to apply, semi-permeable for breathability
- SpeedConnect™ Tubing — low-profile flange, pliable and kink-resistant tubing

White foam dressing — hydrophilic, high density, smaller pore size, flexible wet or dry

EtO sterilized





Ordering information

Item	Size	Cat. No.	Qty.
Black Foam Dressing Kit	Small, 10x8x3cm	47-1702	10 ea/cs
		47-1702/15	15 ea/cs
	Medium, 20x12.5x3cm	47-1701	10 ea/cs
		47-1701/15	15 ea/cs
	Large, 25x15x3cm	47-1700	10 ea/cs
		47-1700/15	15 ea/cs
	XLarge, 58.5x33x3cm	47-1703	10 ea/cs
White Foam Dressing		47-1751	10 ea/cs

CardinalHealth
Essential to care™

Canisters

Features:

- Robust, impact resistant design
- Filters to help prevent leaks and control odor
- Gel packs to solidify exudate
- Easy to install and click into place
- Integrated volume graduation
- Fits discreetly in the device



Ordering information

Item	Description	Cat. No.	Qty.
300cc canister	with gel	47-4000	10 ea/cs
500cc canister	with gel	47-4500	10 ea/cs

Cardinal Health™ NPWT accessory options

Features


Customize your wound care for exactly what your patient needs.

- Bridging Kit
- Y-Connector
- Irrigation tubing with SpeedConnect™ tubing
- Irrigation delivery kit
- SensiSkin™ drape for sensitive skin
- IV Pole adaptors
- Carry cases
- AC power adaptors











Ordering information

Item	Description	Cat. No.	Qty.
Polyurethane Drape		47-7000	10 ea/cs
SensiSkin™ Drape		47-7100	10 ea/cs
Bridging Kit		47-1704	5 ea/cs
Y-Connector		47-2500	10 ea/cs
Irrigation tubing with SpeedConnect™ tubing		47-6000	5 ea/cs
Irrigation delivery kit		47-6500	5 ea/cs
PRO Family	Carrying Bag	47-9600	1 ea/cs
	IV Pole Holder	47-5600	1 ea/cs
	AC Power Supply	47-9100	1 ea/cs

For more information, contact your Cardinal Health Canada Sales representative, or email woundcare@cardinalhealth.ca

© 2015 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, ESSENTIAL TO CARE, SENSISKIN, SPEEDCONNECT, SVEDMAN, SVED and WOUND CARE. ANYWHERE. are trademarks or registered trademarks of Cardinal Health. All other marks are the property of their respective owners. Lit. No. 2MC15-29403 (10/2015). PRC approved (10/15).