

Seeing is Believing with IRIS Technology

Feeding Tube Placement — No Longer a Blind Procedure

Kangaroo™ Feeding Tube with IRIS Technology features a small bore feeding tube with an integrated camera to assist clinicians during feeding tube placement through visualization of key anatomical markers.



CardinalHealth™



Kangaroo™ Console



Kangaroo™ Feeding Tube with IRIS Technology



Kangaroo™ Interface Cable



Kangaroo™ Insufflation Device

Kangaroo™ Console

- 10.4 in. (26.4 cm) touch screen
- Still-shot capture and exportability of saved images
- In-use battery life up to four hours
- Multilingual

Kangaroo™ Feeding Tube with IRIS Technology

- Integrated 3 mm camera at the distal tip of the tube
- 8 Fr, 10 Fr, and 12 Fr sizes available
- ENFit™* Connection

Kangaroo™ Interface Cable

- Still-shot capture button
- 5 ft. (1.5 m) length for easy maneuverability
- Reusable for up to 5000 connections

Kangaroo™ Insufflation Device

- Manual insufflation device to deliver air as needed during placement
- Delivers approximately 30 mL of air per squeeze
- Built-in pressure release valve
- ENFit™* Connection

Ordering Information

Code Number	Description	Case Qty
Kangaroo™ Console		
Z386100	Console with power adapter, battery, interface cable, and pole clamp	1 ea/ca
Kangaroo™ Feeding Tube with IRIS Technology (non-sterile)		
Z460836E	8 Fr x 36 in. (91 cm)	5 ea/ca
Z460843E	8 Fr x 43 in. (109 cm)	5 ea/ca
Z460855E	8 Fr x 55 in. (140 cm)	5 ea/ca
Z461036E	10 Fr x 36 in. (91 cm)	5 ea/ca
Z461043E	10 Fr x 43 in. (109 cm)	5 ea/ca
Z461055E	10 Fr x 55 in. (140 cm)	5 ea/ca
Z461236E	12 Fr x 36 in. (91 cm)	5 ea/ca
Z461243E	12 Fr x 43 in. (109 cm)	5 ea/ca
Z461255E	12 Fr x 55 in. (140 cm)	5 ea/ca
Kangaroo™ Accessories		
Z388103	Interface Cable	1 ea/ca
Z388104	Pole Clamp	1 ea/ca
Z388105E	Insufflation Device	5 ea/ca
Z388102	Replacement Battery	1 ea/ca
Z388101	Global Power Adapter	1 ea/ca

For more information, contact your local sales representative or customer service at 1.888.291.5033.

© 2018 Cardinal Health Canada. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO and Kangaroo are trademarks of Cardinal Health and may be registered in the US and/or in other countries.
Lit. No. H9812, CA-PMR-0158, 18-200 (12/2018)

ENFit™* is a federally registered trademark in multiple jurisdictions throughout the world. Used with permission.

Cardinal Health Canada
1000 Tesma Way
Vaughan, ON L4K5R8
1.888.291.5033
shop.cardinalhealth.ca

CardinalHealth™